

JAY[®]

CONFIGUREFIT™

≡ Tailored seating solutions for the user requiring a custom CONFIGURED product with the ultimate personal FIT.

- ≡ Alignment to Reimbursement:
- ≡ Ease of Ordering:
- ≡ Improved and Expanded Product Line:
- ≡ Single Location Assembly of Seating with Mobility Base:



Ease of Ordering

The ordering process addresses positioning needs in clinical order by region. Each region is divided into zones representing specific positioning support. The zones are followed by configurable item numbers, each with the appropriate HCPCS code. The clinical need, product required and billing indication are all aligned.

Region	Zone	Part Number	Description	HCPCS Code	
Base		Zone 1	CS-01-SEAT	Custom Configured Seat	E2609
		Zone 2	CS-02-PELVICSPU CS-02-PELVICSP CS-02-PELVICDPP	Single Pull Anterior Pelvic Support Unpadded Single Pull Anterior Pelvic Support Padded Dual Pull Anterior Pelvic Support Padded	E0978 E0978 E0978
		Zone 3	CS-03-HIP CS-03-HIP_REM CS-03-MODU	Lateral Pelvic/Thigh Supports with Any Type Hardware Upgrade Lateral Pelvic/Thigh Supports to Removable Hrdw Upgrade Lateral Pelvic Supports to Modular Hardware Lateral Hip Support Modifications	E0956 E1028 K0108 K0108
		Zone 4	CS-04-ADD CS-04-ADD_REM CS-04-MODU CS-04-ADD_FM CS-04-ADD_SA	Lateral Adductor Supports with Any Fixed Hardware Upgrade Lateral Adductor Supports to Removable or Swing Away Upgrade Lateral Adductor Supports to Modular Hardware Frame Mount Lateral Adductor Supports with Fixed Hardware Upgrade Frame Mount Adductor Supports to Swing Away Lateral Adductor Support Modifications	E0956 E1028 K0108 E0956 E1028 K0108
		Zone 5	CS-05-ABD CS-05-ABD_REM	Medial Thigh Support with Fixed Bracket Medial Thigh Support Removable/Retractable Upgrade Medial Thigh Support Custom Modifications	E0957 E1028 K0108
Lower Extremity		Zone 6		Lower Extremity Supports – See Product Selection Guide for Styles and Sizes	E0951/E0952
Torso		Zone 7	CS-07-BACK	Custom Configured Back	E2617
		Zone 8	CS-08-LAT CS-08-LAT_REM CS-08-LAT_SA CS-08-MODU CS-08-1130 CS-08-CF	Lateral Thoracic Supports with Any Type Hardware Upgrade Lateral Thoracic Supports to Removable Hardware Upgrade Lateral Thoracic Supports to Swing Away Hardware Upgrade Lateral Thoracic Supports to Modular Hardware Upgrade Swing Away Supports to Quick Adjust Hardware Upgrade to Contour-Fit Lateral Pad Lateral Thoracic Support Modifications	E0956 E1028 E1028 K0108 K0108 K0108 K0108
		Zone 9	CS-09-ATS1 CS-09-ATS2 CS-09-ATSBF CS-09-ATSTS	Anterior Thoracic Support Center or Zipper Open Style Anterior Thoracic Support Butterfly Style Anterior Thoracic Support Anterior Trunk Strap Support	E0960 E0960 E0960 E0960
		Zone 10	CS-10-PRO	Protractor Supports	K0108
		Zone 11	CS-11-RET	Retractor Supports	K0108
Upper Extremity		Zone 12	CS-12-UES	Arm Troughs with Hardware	E2209
		Zone 13	CS-13-TRAY CS-13-TRAY_MNT	Upper Extremity Support Wheelchair Tray Tray Mounting Hardware	E0950 K0108
Head		Zone 14	CS-14-HEAD CS-14-REM CS-14-MULTI CS-14-AACC CS-14-LACC CS-14-HEAD-OS CS-14-REM-OS CS-14-3P	Head Support with Fixed Hardware Upgrade to Removable Hardware Upgrade to Multi-Pad System Anterior Head Support Accessories Lateral Head Support Accessories Head Support with Fixed Hardware Upgrade to Removable Hardware 3 Panel Head Support Accessories	E0955 E1028 K0108 K0108 K0108 E0955 E1028 K0108

Improved and Expanded Product Line

		JAY ConfigureFit™ Description	Benefit
SEATS	GENERAL	The order process has changed so that any style of seat can be ordered with any style of attaching hardware. There is only one master part number for seats, this part is billed as E2609 and is ordered separately from the back so each back can also be combined with any seat and be billed separately as an E2617 HCPCS Code. Style sheets have been introduced to control the upholstery style. The default style, Style A, produces a finish like AES. To obtain a style like RehabilTech select Style R.	Any seat can be ordered with any hardware. Clearer billing process. No compromise on style of upholstery.
	T NUT PATTERN AND EDGE TREATMENT	Seats have adopted the AES t-nut pattern and will default to Style A (AES upholstery and wood treatment) Style R will allow for wood and upholstery treatments to follow RehabilTech style. This should be indicated on the order form.	A t-nut pattern that is consistent.
	STANDARD FOAM	Standard foam is 1.5" high resiliency foam with override allowed on both the order form and through profile.	A foam standard that is consistent.
	FABRICS AND COLORS	ConfigureFit adopts the AES color pallet with the addition of Black Vinyl (V014) from RehabilTech and 3DX Spacer Fabric.	More fabrics.
	CONTOURED SEATS AND ANTI-THRUST SEATS	ConfigureFit adopts the AES style with Quick Fit Adjustable Contour and Quick Fit Adjustable Anti-Thrust seats.	Maintains Adjustability.
	MODIFICATIONS	The available modifications you have come to expect from custom seating are all available on the ConfigureFit product line as well. They are a combination of modifications from both the AES and RehabilTech product lines.	Expanded Product Line with even more modifications.
	ATTACHING HARDWARE	All Universal Hardware is adopted from the AES product line with new attaching hooks and forks (L brackets) similar to the RehabilTech Style. Add a locking pin to this hardware for transit. ConfigureFit has all new drop hooks with a new angle adjustment lock-out feature and a longer slot. Swivel tabs and stops come from the AES line. ShurLocks, EZ Mount, snap on hardware etc, are all still available. A new snap on clamp similar to Freedom Designs will be offered as well.	Universal hardware eliminates need for a specific clamp to match the chair's tubing diameter (benefit to RehabilTech users) Simpler transit hardware has more consistency, and simpler drop hook hardware (only one style whether angle adjustable or not) simplifies ordering.
	ADJUSTMENT RAILS (TRACKS)	Adjustment rails (or tracks) are now ordered as part of the seat. Single and double styles are available with an option to allow the tracks to cover the full depth of the seat.	Improves margin with cost rolled into E2609, adds Double Slot Rails (Tracks) for AES users.
PELVIC BELTS	STYLES	AES ComforFit products are still available with many of the same part numbers you are used to. Some upper level part numbers and pricing has changed to improve margin on select products. Some select specialty items are retained from the RehabilTech product line.	Maintains a familiar product offering.
	SIZES	No Changes.	
LATERAL PELVIC SUPPORT	GENERAL	Each Pelvic Support is organized as pad and fixed hardware to align to the E0956 HCPCS code. The item can be upgraded to a removable bracket to align to the E1028 HCPCS code. The brackets can be also upgraded to modular style to align to the K0108 HCPCS code.	Improved alignment to HCPCS codes.
	PADS	Pads are offered in standard sizes (for pelvic support only) or as full length or in custom sizes. An up charge must be selected to capture the charge for pads that are not full length or standard sized. Contoured full length pads are no charge. The pads consist of a standard plywood base and 1/2" of high resiliency foam. Pads are ordered as separate left and right with a difference in sizing from left to right allowed.	Left and right selection allows for one HCPCS code per item, rather than billing a quantity of two on the HCPCS item. No up charge for contoured pads.
	FIXED BRACKETS	Fixed brackets are available as standard, heavy duty and 15 degree angled which can be used to taper the distance between the pads from front to rear. All brackets are styled like the AES bracket. Each is available in four profiles (or heights) to accommodate different seat thickness. They are also available in four offsets.	Expanded product line.
	REMOVABLE BRACKETS	Removable Brackets are available as standard removable like the AES style, Rail mount is available for mounting on mobility bases with tracking along the side frame, as well as 20 degree adjustable, (can be used to taper the distance between the pads from front to rear) and Large Target Button (RehabilTech). Each is available in four profiles (or heights) to accommodate different seat thickness. They are also available in four offsets.	Expanded product line set up to align to E1028 HCPCS code.
	MODULAR BRACKETS	Modular and 20 degree adjustable Modular brackets are available from the AES product line. Each is available in four profiles (or heights) to accommodate different seat thickness. They are also available in four offsets.	Aligned to K0108 HCPCS Code.
	MODIFICATIONS	The available modifications you have come to expect from custom seating are all available on the ConfigureFit product line. They are a combination of modifications from both the AES and RehabilTech product lines.	Expanded product line with even more modifications.
LATERAL ADDUCTOR SUPPORT	GENERAL	Each adductor support is organized as pad and fixed hardware to align to the E0956 HCPCS code. The item can be upgraded to a removable bracket to align to the E1028 HCPCS code. The brackets can be also upgraded to modular to align to the K0108 HCPCS code.	Improved alignment to HCPCS codes.
	PADS	Pads are offered in standard sizes for adductor support or in custom sizes. An up charge must be selected to capture the charge for pads that are not standard sized. The pads consist of a standard plywood base and 1/2" of high resiliency foam. Pads are ordered as separate left and right with a difference in sizing from left to right allowed.	Left and right selection allows for one HCPCS code per item, rather than billing a quantity of two on the HCPCS code. No up charge for contoured pads.
	FIXED BRACKETS	Fixed brackets are available as standard, heavy duty and 15 degree angled which can be used to taper the distance between the pads from front to rear. All brackets are styled like the AES brackets. Each is available in four profiles (or heights) to accommodate different seat thickness. They are also available in four offsets.	Expanded product line.
	REMOVABLE BRACKETS	Removable Brackets are available as standard removable like the AES style, Rail mount is available for mounting on mobility bases with tracking along the side frame, as well as 20 degree adjustable, (can be used to taper the distance between the pads from front to rear) and Large Target Button (RehabilTech). Each is available in four profiles (or heights) to accommodate different seat thickness. They are also available in four offsets.	Expanded product line set up to align to E1028 HCPCS code.
	MODULAR BRACKETS	Modular and 20 degree adjustable Modular brackets are available from the AES product line. Each is available in four profiles (or heights) to accommodate different seat thickness. They are also available in four offsets.	Aligned to K0108 HCPCS Code.
	FRAME MOUNT	Frame mount adductors are also available. These consist of a round pad and a fixed (or if upgraded) swing away bracket. These are ordered as an alternative to adductors attached to the seat.	
	MODIFICATIONS	The available modifications you have come to expect from custom seating are all available on the ConfigureFit product line. They are a combination of modifications from both the AES and RehabilTech product lines.	Expanded product line with even more modifications.
MEDIAL THIGH SUPPORT	GENERAL	Each medial thigh support is organized as a pad and fixed hardware to align to the E0957 HCPCS code. The item can be upgraded to a removable bracket to align to the E1028 HCPCS code.	Improved alignment to HCPCS codes.
	PADS	Pads are offered in both oval and wedge designs. Pads come in standard sizes for medial thigh support or in custom sizes. An up charge must be selected to capture the charge for pads that are not standard sized. The pad construction is based on AES design.	Expanded Product Line.
	FIXED BRACKETS	Fixed brackets are available as standard.	Aligned to E0957 code.
	REMOVABLE BRACKETS	Flip Down Brackets are available in many styles from both the AES and RehabilTech product lines. This includes the AES Pull Away, Hide Away and Swing Down styles as well as the RehabilTech Push Button. The RehabilTech push button is available as standard and narrow profile with the removable feature for both.	Expanded product line set up to align to E1028 HCPCS code.
	MODIFICATIONS	The available modifications you have come to expect from custom seating are all available on the ConfigureFit product line. They are a combination of modifications from both the AES and RehabilTech product lines.	Expanded product line with even more modifications.

≡ How will the new ConfigureFit product differ from AES and RehabiliTech?

As Compared to AES
When ordering Style A there will be no difference with the AES Product Line. There are differences with attaching hardware detailed below.
No change.
No change.
Same colors and fabrics with RehabiliTech Black Vinyl added and 3DX Spacer Fabric added.
No change.
No change.
Universal Hardware is the same, however new attaching hooks and forks are utilized that will be similar to the RehabiliTech style. All old transit hardware has been discontinued. The new hooks and forks require the addition of a locking pin for transit. The drop hooks are all new, with the old style discontinued. Hardware is now bundled into kits for mounting at the front or rear of the seat. These are selected individually for a mix and match functionality.
Double slot tracks are now available. No other changes.
Continued ComfortFit offering with the addition of some RehabiliTech specialty items.
No change.
All pelvic support items are similar in concept to Basis packages for alignment to HCPCS codes.
Construction is unchanged. Some standard sizing changes.
Adds heavy duty steel to the offering.
Add RehabiliTech Removable.
No change.
No change.
All adductor support items are similar to Basis packages for alignment to HCPCS codes.
Construction is unchanged. Some standard sizing changes.
Adds heavy duty steel to the offering.
Add RehabiliTech Removable.
No change.
No change.
No change.
All medial thigh support items are similar to Basis packages for alignment to HCPCS codes.
Construction is unchanged. .
No change.
Adds RehabiliTech push button and removable.
No change.

As Compared to RehabiliTech
When ordering Style R there is little difference from the RehabiliTech Product Line with minor T-nut pattern changes. There are differences with attaching hardware detailed below.
Small Change. A new ConfigureFit Seat will not mate to an old RehabiliTech back if seat to back brackets were used. The t-nut alignment is not compatible. Must order both seat and back as new with new seat to back bracket.
1/2" thicker than RehabiliTech standard.
Adopts all AES colors and fabrics plus Black Vinyl from RehabiliTech (V014) and 3DX Spacer Fabric.
To obtain old style contoured seat, must now order Quick Fit Adjustable Contour. It has some differences but maintains function. Anti-thrust seat adds additional benefit of Quick Fit adjustability when this is specified.
No change.
The clamp portion of knob release hardware is derived from the AES style universal hardware. Forks are essentially the same. Drop hooks are all new with angle adjustability. These can be locked to a 90 degree orientation with a built in lock-out feature. Clamp on hook stops are adopted from the AES line and eliminates the spacer found on the carriage bolt of RehabiliTech clamps. Hardware is now bundled into kits for mounting at the front or rear of the seat. These are selected individually for a mix and match functionality.
RehabiliTech style single slot tracks are discontinued, replaced, instead, by the AES style single slot track. For extra security order the double slotted track.
Continued ComfortFit offering with the addition of some RehabiliTech specialty items.
No change.
All pelvic support items are similar in style to 18E0956 and 18E1028X packages for alignment to HCPCS codes.
Construction is unchanged. An introduction of standard sizing, no charge for contoured pads.
The brackets are styled like AES brackets. The offering is expanded with 15 degree brackets and heavy duty has been lightened considerably.
Adds AES Removable, Rail Mount and 20 degree adjustable.
Expanded product line.
No Change
All adductor support items are similar to 18E0956 and 18E1028X packages for alignment to HCPCS codes
Construction is unchanged. An introduction of standard sizing.
The brackets are styled like AES brackets. The offering is expanded with 15 degree brackets and heavy duty has been lightened considerably.
Adds AES Removable, Rail Mount and 20 degree adjustable.
Expanded product line.
No change.
No change.
All medial thigh support items are similar to 24E0957 and 24E1028 packages for alignment to HCPCS codes.
Construction is based on AES. An introduction of standard sizing.
No change.
Adds AES Hide Away, Pull away and Swing Down. RehabiliTech Swing Down is discontinued.
No change.

Improved and Expanded Product Line (Continued)

		JAY ConfigureFit™ Description	Benefit
LOWER EXTREMITY SUPPORTS	SHOE HOLDERS	All shoe holders available from both AES and RehabilTech product lines. The items are set up as left and right ordering.	All shoe holder products still available. Left and right ordering allows easy mix and match sizing.
	FOOT AND ANKLE POSITIONERS	All foot and ankle positioners from the ComforFit line up are available, coded as K0108.	Better HCPCS Code alignment.
	FOOT BOXES	Foot Boxes adopted from the AES product line. Some RehabilTech items added to this product category. Many options are available.	Expanded product line.
	MODIFICATIONS	The available modifications you have come to expect from custom seating are all available on the ConfigureFit product line. They are a combination of modifications from both the AES and RehabilTech product lines.	Expanded product line with even more modifications.
BACKS	GENERAL	The order process has changed so that any style of back can be ordered with any style of attaching hardware. There is only one master part number for backs, which is billed as E2617 and is ordered separately from the seat so each back can also be combined with any seat and then billed separately as an E2609. Style sheets have been introduced to control the upholstery style. The default style, Style A, produces a finish like AES. To obtain a style like RehabilTech select Style R.	Any back can be ordered with any hardware. Clearer billing process. No compromise on style of upholstery.
	T NUT PATTERN AND EDGE TREATMENT	Backs have adopted the AES t-nut pattern and will default to Style A (AES upholstery and wood treatment) Style R will allow for wood and upholstery treatments to follow RehabilTech style. This should be indicated on the order form.	A T-nut pattern that is consistent.
	STANDARD FOAM	Standard foam is 1" high resiliency foam with an override allowed on both the order form and through profile.	A foam standard that is consistent.
	FABRICS AND COLORS	ConfigureFit adopts the AES color pallet with the addition of Black Vinyl (V014) from RehabilTech and 3DX Spacer Fabric.	More fabrics.
	OTHER	Bi angular backs are available with many bracket types from both the AES and RehabilTech product lines. Contoured backs and Foam Form backs are available and are derived from the AES product line.	Expanded product line.
	MODIFICATIONS	The available modifications you have come to expect from custom seating are all available on the ConfigureFit product line. They are a combination of modifications from both the AES and RehabilTech product lines.	Expanded product line with even more modifications.
	ATTACHING HARDWARE	Universal Hardware has been adopted from the AES product line with new attaching hooks and forks similar to the RehabilTech style. Add a locking pin for transit. New adjustable angle forks (used with universal hardware) are available for curved backs. ConfigureFit has all new drop hooks with a new angle adjustment lock-out feature and a longer slot. Swivel tabs and stops come from the AES line. ShurLocks, EZ Mount, snap on hardware etc, are all still available.	Universal hardware eliminates need for a specific clamp to match tubing diameter (Benefit to RehabilTech users) Simpler transit hardware with more consistency, simpler drop hook hardware (only one style whether angle adjustable or not) Simplified hardware for curved backs which eliminates many bracket angles in favor of an adjustable angle bracket.
	ADJUSTMENT RAILS (TRACKS)	Adjustment rails (or tracks) are now ordered as part of the back. Single and double styles are available with an option to allow the tracks to cover the full height of the back.	Improves margin with cost rolled into E2617, adds Double Slot Rails (tracks) for AES users.
LATERAL THORACIC SUPPORT	GENERAL	Each thoracic support is organized as a pad and fixed hardware to align to the E0956 HCPCS code. The item can be upgraded to a removable bracket to align to the E1028 HCPCS code. The brackets can be also upgraded to modular to align to the K0108 HCPCS code.	Improved alignment to HCPCS codes.
	PADS	Pads are offered in standard sizes or in custom sizes. An up charge must be selected to capture the charge for pads that are not standard sized. Pads can be curved, or flat, D Shaped or tapered. Base materials can be wood, aluminum or plastic with choice of closed cell foam (1/4") or high resiliency foam (1/2"). Pads are ordered as separate left and right with a difference in sizing from left to right allowed. Pads can also be ordered with an additional left and/or right for three point positioning or a four pad support setup.	Left and right selection allows for one HCPCS code per item, rather than billing a quantity of two on the HCPCS item. Many added options for foam and base materials.
	FIXED BRACKETS	Fixed brackets are available as standard, heavy duty and 15 degree angled, which can be used to taper the distance between the pads from top to bottom. All brackets are styled like the AES bracket. Each is available in four profiles (or heights) to accommodate different seat thickness. They are also available in four offsets. Adjustable profile (height) brackets are derived from the AES product line and can be ordered for an additional charge.	Expanded product line.
	SWING AWAY BRACKETS	Standard swing away brackets are available from the AES product line. 20 degree adjustable are still available. A new angle adjustable swing away is now offered that can replace the AES style or can be ordered as a proximal style to replace the old RehabilTech version. Each is available in different profiles (or heights) to accommodate different back thickness. They are also available in four offsets.	New angle adjustable swing away lateral bracket that can be set up standard (like AES modular) or as proximal (like RehabilTech). Expanded product line set up to align to E1028 HCPCS code.
	MODULAR BRACKETS	Modular and 20 degree adjustable modular brackets are available from the AES product line. Each is available in four profiles (or heights) to accommodate different seat thickness. They are also available in four offsets.	Aligned to K0108 HCPCS Code.
	REAR QUICK ADJUST (SUMMER/WINTER)	Rear Quick Adjust is available as well as Lever style from the RehabilTech product line.	Expanded product line.
	MODIFICATIONS	The available modifications you have come to expect from custom seating are all available on the ConfigureFit product line. They are a combination of modifications from both the AES and RehabilTech product lines.	Expanded product line with even more modifications.
ANTERIOR THORACIC	STYLES	AES ComforFit products are still available in the part numbers you are used to. Some pricing has changed to improve margin on select products. Some select specialty items are retained from RehabilTech product line.	Maintains a familiar product offering.
	SIZES	No changes.	
LAP TRAYS	SIZES	AES sizes have been adopted. Sizes are only slightly different from RehabilTech sizes. Custom sizes are still available. All materials are still available as well.	Maintains a familiar product offering.
	HARDWARE	Hardware is a combination of RehabilTech and AES hardware.	Expanded product line.
HEAD SUPPORTS	GENERAL	Each head support is organized as a pad and fixed hardware to align to the E0955 HCPCS code. The item can be upgraded to a removable or flip down bracket to align to the E1028 HCPCS code. Pads are single pads only. Upgrading to multi-pads requires additional pick to align to the K0108 HCPCS code.	Improved alignment to HCPCS codes.
	PADS	Pads are offered from the Whitmyer product line and some from the AES and RehabilTech product lines. All pads are sold first (under the E0955 HCPCS codes) as single pads. Adding additional pads requires selections that will fall under the K0108 HCPCS code.	Better Alignment to HCPCS Codes, familiar products are all still available.
	BRACKETS	All brackets from the Whitmyer line continue to be available. Brackets are set to the E1028 HCPCS code.	Pricing better balanced between E0955 and E1028 HCPCS Codes.
	OLD STYLE REHABILITECH AND AES HEAD SUPPORTS	ConfigureFit makes some legacy AES and RehabilTech style brackets and pads available under this special section. These pads (such as the three panel) and brackets (such as the anterior/posterior adjustable bracket) are a hybrid of RehabilTech and AES items.	Expanded Product Line - maintains legacy products.
	ACCESSORIES	Accessories are organized under the K0108 HCPCS Code.	Better alignment to HCPCS codes. Familiar products are all still available.
	MODIFICATIONS	The available modifications you have come to expect from custom seating are all available on the ConfigureFit product line. They are a combination of modifications from both the AES and RehabilTech product lines.	Expanded product line with even more modifications.

≡ How will the new ConfigureFit product differ from AES and RehabiliTech?

As Compared to AES	As Compared to RehabiliTech
No change.	No change.
No change.	Add Foot and Ankle Positioners.
No change.	Adopt AES style foot boxes.
No change.	No change.
When ordering Style A there will be no difference with the AES Product Line. There are differences with attaching hardware detailed below.	When ordering Style R there is little difference from the RehabiliTech Product Line with minor T-nut pattern changes. There are differences with attaching hardware detailed below.
No change.	Small change. A new ConfigureFit seat will not mate to an old RehabiliTech back if seat to back brackets are used. The T-nut alignment is not compatible. Must order both seat and back as new with new seat to back bracket.
No change.	No change.
Same colors and fabrics with RehabiliTech Black Vinyl added and 3DX Spacer Fabric added.	Adopts all AES colors and fabrics plus Black Vinyl from RehabiliTech (V014) and 3DX Spacer Fabric.
More options added.	More options added.
No change.	No change.
Universal Hardware is the same, however new attaching hooks and forks are utilized that will be similar to the RehabiliTech style. All old transit hardware has been discontinued. All old back hooks for curved wood backs have been discontinued and replaced by a new angle adjustable hook. The new hooks and forks require the addition of a locking pin for transit. The drop hooks are all new, and the old style has been discontinued. Hardware is now bundled into kits for mounting at the front or rear of the seat. These are selected individually for a mix and match functionality.	The clamp portion of knob release hardware is derived from AES style universal hardware. Hooks are essentially the same. Drop hooks are all new with angle adjustability. These can be locked to a 90 degree orientation with a built in lock-out feature. Clamp on hook stops are adopted from the AES line and eliminates the spacer found on the carriage bolt of RehabiliTech clamps. Hardware is now bundled into kits for mounting at the front or rear of the seat. These are selected individually for a mix and match functionality.
Double slot tracks are now available. No other changes.	RehabiliTech style single slot tracks are discontinued, instead using the AES style single slot track. For extra security order the double slotted track.
All thoracic support items are similar in function to Basis packages for alignment to HCPCS codes.	All thoracic support items are similar in function to 20E0956 and 20E1028X packages for alignment to HCPCS codes.
Construction is unchanged if ordering high resiliency foam. Some standard sizing changes.	Construction is unchanged if ordering closed cell foam. An introduction of standard sizing, with more standard base and foam options.
Adds heavy duty steel to the offering.	The brackets are styled like AES brackets. The offering is expanded with 15 degree brackets and heavy duty has been lightened considerably. Adjustable height brackets are discontinued with the AES bracket utilized instead.
Add new angle adjustable bracket with proximal option.	Discontinue RehabiliTech swing away lateral bracket and replace with a new angle adjustable proximal design. Addition of AES style brackets.
No change.	Expanded product line.
Front Quick Adjust discontinued. Lever style from RehabiliTech now available.	Add Rear Quick Adjust for select swing away and modular laterals.
No change.	No change.
Continued ComforFit offering with addition of some RehabiliTech specialty items.	Continued ComforFit offering with addition of some RehabiliTech specialty items.
No change.	No change.
Addition of 1/4" Polycarbonate trays	Utilizes AES standard sizing.
No change.	Expanded product line.
All head support items are similar in concept to Basis packages for alignment to HCPCS codes.	All head support items are similar in concept to 26E0955 and 26E1028X packages for alignment to HCPCS codes.
All products still available.	All products still available.
No change.	No change.
Add new anterior/posterior bracket.	Discontinued RehabiliTech straight and offset brackets, replaced by AES style. Reinroduce the anterior/posterior style bracket under a new design.
No change.	No change.
No change.	No change.

JAY ConfigureFit

NEW OPTIONS

Expanded Swing-away Lateral Thoracic Bracket Selection:

- Now optionally available with angle adjustability in both standard and proximal configurations.

Transit Option:

- WC20 Transit Tested



SINGLE LOCATION ASSEMBLY OF SEATING WITH MOBILITY BASE

Increased Efficiency Throughout the Process:

- Improves and ensures compatibility between the mobility base and seating device before delivery
- No shipping parts from factory to factory
- Improved delivery time

SPECIFICATIONS

JAY®

JAY ConfigureFit	
System Weight	Dependent upon system size and componentry
Product Size	Seat Width - Minimum 10", all other dimensions variable
Available Foam Materials	Viscoelastic, Sunmate and High Resiliency Foams
Available Cover Materials	3DX, Vinyls, Nylons, Neoprenes, Lycras and other Four-way Stretch Fabrics
Fabric Colors	Wide selection - refer to fabric/color samples
Product Configuration	Seat Only, Back Only, Seat and Back - all with available accessories, or Accessories Only
User Weight Capacity	Non-transit: 250 lbs. Non-transit HD: 300 lbs. Transit: Minimum 48 lbs. - Maximum 200 lbs.
Warranty	Sunrise Limited Warranty