

Maxi-Drain Transfer Bench

- Provides assistance to those who have difficulty stepping over bathtub walls.
- Anodized aluminum frame is lightweight, durable and rust-resistant.
- Drainage holes in seat minimize water build-up.
- Side support arm aids in lateral transfer.
- Seat height is adjustable in ½" increments.
- Reversible back attachment to accommodate any bathroom.
- Attractive platinum grey color.
- Comes assembled with tool-free arm and back attachments.
- 350 lb maximum weight capacity.
- Designed to meet the requirements of HCPCS code: E0247.
- Lifetime limited warranty.

7927A 2/cs



Item	Overall Width	Overall Depth	Seat Width	Seat Depth	Minimum Seat Height	Maximum Seat Height	Clearance Beneath Frame	Height of Back Support Above Seat	Weight Capacity
7927A	29 ¼"	20"	26 ½"	16"	18 ¼"	22 ¼"	15 ¼" - 19 ¼"	15"	350 lb



Graham-Field Health Products

2935 Northeast Parkway
Atlanta, Georgia 30360

2007© Graham-Field Health Products

For Placing Orders:

1-800-347-5678

Fax: 1-800-726-0601

www.grahamfield.com

GF0800039RevA08