

Padded Transfer Bench

- Provides assistance to those who have difficulty stepping over bathtub walls.
- Comfortable, cushioned seat and backrest are water-tight and easy to clean.
- Anodized aluminum frame is lightweight, durable and rust-resistant.
- Side support arm aids in lateral transfer.
- Seat height is adjustable in 1" increments.
- Reversible to accommodate any bathroom.
- Attractive platinum grey color.
- Comes assembled with tool-free arm and back attachments.
- 300 lb maximum weight capacity.
- Designed to meet the requirements of HCPCS code: EO247.
- Lifetime limited warranty.

7956A Commode Transfer Tub Bench 1ea
Comes complete with commode pail and cover.



ITEM	Overall Width	Overall Depth	Seat Width	Seat Depth	Minimum Seat Height	Maximum Seat Height	Clearance Beneath Frame	Height of Back Support Above Seat	Diameter Commode Opening	Weight Capacity
7956A	29 1/4"	16 3/4"	25"	15 3/4"	17 1/2"	22"	14" - 18 1/2"	16"	8"	300 lb



Graham-Field Health Products

2935 Northeast Parkway
Atlanta, Georgia 30360

2005© Graham-Field Health Products

For Placing Orders:

1-800-347-5678

Fax: 1-800-726-0601

www.grahamfield.com