



Bathboard with handle

For ease of transfer and enhanced user stability, this moulded plastic bathboard comes complete with a handle. The bathboard provides a comfortable seating area and has holes to aid drainage. The width adjustable brackets are made from non-slip rubber to give a secure fit.

Bathboard with handle 10160C

Technical specification

Length	690 mm (27")	Product weight	1.5 kg
Width at narrowest point	260 mm (10.25")	Max user weight	113 kg (18 st)



Merlin bathseat

A strong bathseat made from reinforced PVC. Fixes securely with large suction pads which incorporate easy-release tabs. Finished in brilliant white with soft grey trim.

Height 152 mm (6") 10157C
Height 203 mm (8") 10158C
Height 305 mm (12") 10159C

Technical specification

Width	280 mm (11")	Product weight	2.5 kg
Length	457 mm (18")	Max user weight	160 kg (25 st)