




Merlin bathboard

Made of highly durable reinforced PVC the Merlin bathboard is compatible with plastic, steel and cast iron baths and is secured into position with a self locating fixing system. The slats are tapered to encourage rapid drainage.

Length 660 mm (26")	10164C
Length 686 mm (27")	10165C
Length 711 mm (28")	10166C

Technical specification

Width	228 mm (9")	Max user weight	190 kg (30 st)
Product weight	2.2 kg		



Merlin showerboard

A strong showerboard constructed of reinforced PVC for maximum durability. Secured into position with a self-locating fixing system. The slats are tapered to encourage rapid drainage.

Length 660 mm (26")	10167C
Length 686 mm (27")	10169C
Length 711 mm (28")	10168C
Optional handle for Merlin bath and showerboards	10176C

Technical specification

Width	340 mm (13.5")	Max user weight	190 kg (30 st)
Product weight	3.2 kg		