

Grab rails



Merlin grab rails

An extra strong grab rail attractively finished in brilliant white. The ribbed handle is 38mm wide to provide the user with increased comfort and grip.

Length 305 mm (12")	10274
Length 457 mm (18")	10276
Length 609 mm (24")	10278

Technical specification

Product weight	375 - 535 g	Diameter	38 mm (1.5")
Depth	76 mm (3")		



Folding support arm

A double height rail for easy rising. Folds up out of the way when not in use. Rail diameter 32 mm (1.25").

Folding support arm 760 mm (30")	10328C
Folding support arm 550 mm (21.5")	10327C

Technical specification

Overall width	125 mm (5")	Max user weight	127 kg (20 st)
Overall height	305 mm (12")		