


B Plastic Grab Rail 305mm
Order code - G26

Length between centres - 305mm (12")
Depth from rail to wall - 80mm (3 1/4")


Max User Wt
157Kg 25st



C Plastic Grab Rail 457mm
Order code - G27

Length between centres - 457mm (18")
Depth from rail to wall - 80mm (3 1/4")

Max User Wt
125Kg 20st




D Plastic Grab Rail 610mm
Order code - G28

Length between centres - 610mm (24")
Depth from rail to wall - 80mm (3 1/4")

Max User Wt
100Kg 16st

Moulded Grab Rail Features

- Ergonomically designed handle with flexible soft-feel undergrip
- Oval plastic extrusion with shaped brackets and concealed screw covers
- Plastic construction eliminates need for earthing
- Maximum width of profile 40mm



Inside grip



Moulded Bracket

E Dropdown Rail
Order code - G29

From wall to end of rail:
Length of rail - 700mm (27 1/2")
Mounting bracket - 260mm (10 1/4") x 100mm (4")
Maximum distance from wall when upright - 145mm (7 3/4")

Max User Wt
125Kg 20st



Warning - Grab Rails

These products should only be fixed to walls of sound construction and careful consideration should be given when choosing fixings.

F Dropdown Rail with Adjustable Leg
Order code - G30

As G29 (Page Listing 'E') but:
Maximum distance from wall when upright - 100mm (4")
Adjustable leg - 690 (27") min & 780mm (31") max

Max User Wt
190Kg 30st



Drop Down Grab Rail Features

A high quality, all steel grab rail with several innovative features:

- Main rail is manufactured from oval section white nylon modified epoxy powder coated steel tube with grey soft feel grey end cap
- 5 hole fixing

- Secondary under-slung supporting arm gives ultimate strength to this stylish extremely robust rail
- Wall mounted bracket is cast in aluminium alloy
- Parks easily in upright position with positive location

