

Platinum Collection Bath Seats

- Anodized aluminum frame is lightweight, durable and rust-resistant
- Seat height is adjustable in 1" increments
- Attractive platinum grey color
- Available with or without backrest
- Comes assembled with tool-free back attachment
- Includes laminated color operating instructions
- Designed to meet the requirements of HCPCS code: E0245
- 350 lb weight capacity
- Lifetime limited warranty

7921A-4 with backrest 4/cs
7931A without backrest 4/cs

Specifications

Overall Width 20"
 Overall Depth 17"
 Seat Width 20"
 Seat Depth 12"
 Minimum Seat Height 14"
 Maximum Seat Height 22"
 Weight Capacity 350 lb

