

Platinum Collection Deluxe Padded Bath Seats

- Comfortable, cushioned seat and backrest are water-tight and easy to clean.
- Anodized aluminum frame is lightweight, durable and rust-resistant.
- Seat height is adjustable in 1" increments.
- Attractive platinum grey color.
- Available with or without support arms.
- Comes assembled with tool-free arm and back attachments.
- 300 lb maximum weight capacity.
- Designed to meet the requirements of HCPCS code: E0245.
- Lifetime limited warranty.

7944A Without Support Arms 1 ea



ITEM	Overall Width	Overall Depth	Seat Width	Seat Depth	Minimum Seat Height	Maximum Seat Height	Height of Back Support Above Seat	Width Between Arms	Weight Capacity
7944A	20 1/2"	21"	15 3/4"	16"	17 1/2"	22 1/2"	16"	NA	300 lb



Graham-Field Health Products

2935 Northeast Parkway
Atlanta, Georgia 30360

2005© Graham-Field Health Products

For Placing Orders:

1-800-347-5678

Fax: 1-800-726-0601

www.grahamfield.com