

Platinum Collection Deluxe Padded Commode Bath Seats

- Open front seat provides easier access for personal care.
- Comfortable, cushioned seat and backrest are water-tight and easy to clean.
- Anodized aluminum frame is lightweight, durable and rust-resistant.
- Seat height is adjustable in 1" increments.
- Attractive platinum grey color.
- Comes complete with commode pail and cover.
- Available with or without support arms.
- Comes assembled with tool-free arm and back attachments.
- 300 lb maximum weight capacity.
- Designed to meet the requirements of HCPCS code: E0245.
- Lifetime limited warranty.

7947A With Support Arms 1 ea



ITEM	Overall Width	Overall Depth	Seat Width	Seat Depth	Minimum Seat Height	Maximum Seat Height	Height of Back Support Above Seat	Diameter Commode Opening	Width Between Arms	Weight Capacity
7947A	20 1/2"	21"	15 3/4"	16"	17 1/2"	22 1/2"	16"	8"	17 3/4"	300 lb



Graham-Field Health Products

2935 Northeast Parkway
Atlanta, Georgia 30360

2005© Graham-Field Health Products

For Placing Orders:

1-800-347-5678

Fax: 1-800-726-0601

www.grahamfield.com