



Self propelled shower/ commode chair

This durable three-in-one self propelled chair can be used as a shower chair, a commode chair or over a toilet pedestal. The non-rust aluminium frame is height adjustable and comes with detachable, height adjustable footrests. The seat is padded for extra comfort and features a removable section for easier personal cleaning. Handles at the rear of the chair make it easy if assistance is required. Padded armrests flip up and out of the way. The two swivelling front castors are 125mm in diameter. The quick release rear wheels are 600mm in diameter and each features a locking hand brake. Supplied with a commode pan.

Self propelled shower/commode chair	10210C
Spare commode pan with lid	10411C

Technical specification

Product weight	17.5 kg	Seat dimensions	430 x 430 mm (17 x 17")
Maximum depth	1140 mm (45")	Seat height	540 - 640 mm (21.25 x 25")
Overall width	610 mm (24")	Max user weight	120 kg (19 st)



Attendant propelled shower/ commode chair

Similar to the self propelled shower/commode chair but supplied with four swivelling castors, each with an integral brake.

Attendant propelled shower/commode chair	10211C
Spare commode pan with lid	10411C

For added convenience the rear wheels are interchangeable between the 125mm and 600mm versions. To convert the attendant propelled to the self propelled version you will need to purchase a conversion kit, consisting of the quick release wheels, axles and hand brakes. To convert the self propelled to the attendant propelled version you will need to purchase a pair of 125mm castors.

Conversion kit	WCAS1488
125 mm castors, pair	WCAS1487

Technical specification

Product weight	14 kg	Seat dimensions	430 x 430 mm (17 x 17")
Maximum depth	1060 mm (41.75")	Seat height	540 - 640 mm (21.25 x 25")
Overall width	540 mm (21.25")	Max user weight	120 kg (19 st)