


Shower stretcher

**LI2403.150, LI2403.170 and LI2403.190**

- HT-cover material
- hinges to the wall (folds up or folds down)
- hinged leg supports
- including wall clamp to secure stretcher when folded up against the wall
- available with adjustable backrest by means of a gas spring (LI2404.150, LI2404.170 and LI2404.190)
- available with pvc water receiver and flexible drain hose (LI2407.150, LI2407.170 and LI2407.190)
- available with adjustable backrest by means of a gas spring and with pvc water receiver and flexible drain hose (LI2408.150, LI2408.170 and LI2408.190)

	A	B
LI2403.150	1400	185
LI2403.170	1090	285
LI2403.190	1800	385




Measurement drawings of LI2404, LI2407 and LI2408 available upon request or see www.linido.com


Height adjustable shower trolley

**LI2514.150, LI2514.170 and LI2514.190**

- HT-cover material
- four stainless steel double braked castors Ø 125 mm
- adjustable backrest by means of a gas spring
- pvc water receiver with flexible hose of 1000 mm
- height adjustable via hydraulic pump
- two hinged side fences
- only available in the mentioned dimensions
- available without backrest (LI2513.150, LI2513.170 and LI2513.190)




	A	B
E514150	1505	465
E514170	1005	615
E514190	1905	615



Removable & hinged side fence

A: Removable side fence LI2490.125

- suitable for all Linido shower stretchers and shower trolleys
- including plastic mounting clamps

**B: Hinged side fence LI2590.125**

- suitable for all Linido shower trolleys
- including plastic mounting clamps

