

# Shower chairs



## Wall mounted shower seat

Corrosion resistant fold up shower seat. The comfortable seat has holes to aid drainage. Available with or without legs.

**NOTE:** Wall mounted shower seats should be fitted by a competent person using screws and rawplugs appropriate for the wall. For this reason we do not include screws or rawplugs with our seats.

Wall mounted shower seat	10242C
Wall mounted shower seat with legs	10243C

### Technical specification

Product weight	3.5 kg	Seat width	495 cm (19.5")
Product weight with legs	5.0 kg	Seat height	405 - 560 mm (16-22")
Seat depth	305 cm (12")	Max user weight	100 kg (16 st)



## Folding shower chair

A folding shower chair with arms and back. Height adjustable. Seat has a cut-away for improved personal hygiene and holes to aid drainage. Folds away when not in use. Supplied with non-slip rubber feet.

Folding shower chair	10178C
----------------------	--------

### Technical specification

Product weight	4.7 kg	Seat width	410 mm (16")
Overall depth	415 mm (16.5")	Seat height	460 - 540 mm (18-21.25")
Seat depth	400 mm (15.75")	Max user weight	113 kg (18 st)
Overall width	530 mm (21")		



## Adjustable height round shower stool

A height adjustable round shower stool with moulded plastic seat and lightweight aluminium legs. Supplied with non-slip rubber feet.

Adjustable height round shower stool	10231C
--------------------------------------	--------

### Technical specification

Product weight	1.5 kg	Seat height	400 - 500 mm (15.75-19.75")
Seat depth	380 mm (15")	Max user weight	113 kg (18 st)
Seat diameter	320 mm (12.5")		