

COOPERS®

Shower chairs



Extra heavy duty shower chair

The reinforced frame and extra wide seating surface mean this shower chair can accommodate users of up to 227 kg (36 st). Supplied with non-slip rubber feet.

Extra heavy duty shower chair 10220C

Technical specification

Product weight	6 kg	Seat height	380 - 480 mm (15-19")
Seat depth	380 mm (15")	Max user weight	227 kg (36 st)
Seat width	510 mm (20")		



Shower mat, pack of 10

Made from natural rubber, this shower mat is machine washable at low temperatures on a delicate cycle. White.

Shower mat, pack of 10 10187

Technical specification

Length	560 mm (22")	Width	560 mm (22")
--------	--------------	-------	--------------



Turner cubicle seat

A foldaway cubicle seat with legs to support the weight of the user safely. Plastic coated for anti-corrosion protection. Height adjustable.

NOTE: Wall mounted shower seats should be fitted by a competent person using screws and rawlplugs appropriate for the wall. For this reason we do not include screws or rawlplugs with our seats.

Turner cubicle seat 10240

Technical specification

Product weight	5.5 kg	Seat height	457 - 584 mm (18-21")
Seat width	457 mm (18")	Max user weight	190 kg (30 st)
Seat depth	330 mm (13")		