

Three-In-One Aluminum Commodes

- 3-in-1 design is versatile and economical.
- Aluminum frame is lightweight and easy to transport.
- Contoured arm rests enhance user comfort.
- Sturdy plastic seat snaps-on easily.
- Comes complete with commode pail, cover and splash shield.
- Limited lifetime warranty.
- 4" rear swivel, locking casters.
- Includes laminated color operating instructions.
- Designed to meet the requirements of HCPCS code: E0164.

2215A With Back Bar and Casters 1 ea
2215A-2 With Back Bar and Casters 2/cs



ITEM	Width Between Arms	Overall Width	Seat Width	Seat Depth	Seat Height Range	Weight Capacity
2215A, 2215A-2	16 3/4"	21 1/4"	14"	16 3/4"	21 1/2"	300 lb



Graham-Field Health Products

2935 Northeast Parkway
Atlanta, Georgia 30360

2005© Graham-Field Health Products

For Placing Orders:

1-800-347-5678

Fax: 1-800-726-0601

www.grahamfield.com