

## Locking Raised Toilet Seat With Removable Armrests

- Designed for those who have difficulty sitting down or standing up from the toilet.
- Adds 4 1/2" to toilet seat height, reducing bending effort.
- Front clamping mechanism locks securely to the toilet and prevents shifting of the seat.
- Wide, contoured seating surface is comfortable.
- White, heavy duty plastic construction is strong and durable.
- Push buttons allow for easy removal of armrests.
- 300 lb maximum weight capacity.
- Designed to meet the requirements of HCPCS code: E0244.
- Lifetime limited warranty.

6487RA 2/cs



Item	Width Between Arms	Overall Width	Seat Height	Seat Depth	Weight Capacity
6487RA	18 1/2"	20 3/4"	4 1/2"	16 1/4"	300 lb



**Graham-Field Health Products**

2935 Northeast Parkway  
Atlanta, Georgia 30360

2007© Graham-Field Health Products

For Placing Orders:

**1-800-347-5678**

**Fax: 1-800-726-0601**

**www.grahamfield.com**

GF0700097RevA07