


Silver Collection

Three-In-One Steel Commodes

- Versatile 3-in-1 design can be used as a bedside commode, raised toilet seat or toilet safety frame.
- Welded steel construction provides strength and stability.
- Durable, grey-coated steel frame is attractive and easy to clean.
- Backbar detaches to accommodate use as toilet safety frame.
- Sturdy plastic seat snaps-on easily.
- Seat height readily adjusts in 1 inch increments.
- Comes complete with commode pail, cover and splash shield.
- Includes laminated color operating instructions.
- Designed to meet the requirements of HCPCS code: E0163.
- Limited lifetime warranty.
- Contoured armrests enhance user comfort.
- Seat height readily adjusts from 15 1/2" to 23 1/2", in 1 inch increments.

7103A-4

4/cs



ITEM	Width Between Arms	Overall Width	Seat Width	Seat Depth	Seat Height Range	Weight Capacity	Replacement Tip
7103A-4	17 1/2"	24"	14"	15"	15 1/2" - 23 1/2"	300 lb	9290



Graham-Field Health Products

2935 Northeast Parkway
Atlanta, Georgia 30360

2005© Graham-Field Health Products

For Placing Orders:

1-800-347-5678

Fax: 1-800-726-0601

www.grahamfield.com