


Imperial Collection Three-In-One Commodes
Steel Three-In-One Commode

- Extra wide platform seat is comfortable for long periods of sitting.
- Commode opening is contoured and indented in the front to accommodate the needs of larger users.
- Welded steel construction and extra wide base provide strength and stability.
- Cushioned back bar enhances user comfort.
- Removable backbar facilitates use as a bedside commode, raised toilet seat or toilet safety frame.
- Comes complete with commode pail, cover and splash shield.
- Includes laminated color operating instructions.
- Designed to meet the requirements of HCPCS code: E0168.
- Limited lifetime warranty.
- Seat height readily adjusts from 16 1/2" to 23", in 1 inch increments.
- 400 lb maximum weight capacity.

7446A 1 ea
7446A-2 2/cs

Bariatric



ITEM	Width Between Arms	Overall Width	Seat Width	Seat Depth	Seat Height Range	Weight Capacity	Replacement Tip
7446A, 7446A-2	23 1/2"	26"	23"	18"	16 1/2" - 23"	400 lb	9140



Graham-Field Health Products
2935 Northeast Parkway
Atlanta, Georgia 30360
2005© Graham-Field Health Products

For Placing Orders:
1-800-347-5678 Fax: **1-800-726-0601**
www.grahamfield.com