


# Permapad



The Permapad pressure reducing cushion is designed to be used as an introductory product for a pressure ulcer prevention plan. It is a cost effective, one piece foam cushion with a contoured surface area designed to ensure high levels of comfort, support and pressure reduction for Low to Medium Risk patients.

## Technical Specifications

<b>Product Code</b>	<b>PL1077 - 17" x 17"</b> <b>PL1078 - 20" x 20"</b>
<b>Maximum Weight Limit</b>	17 stone / 108 kgs
<b>Risk Factor</b>	Low to Medium Risk*
<b>Warranty</b>	2 years
<b>Cushion Dimensions</b>	<b>PL1077 - Std. Size: 43cm x 43cm x 10cm depth / 17" x 17" x 4" depth</b> <b>PL1078 - Std. Size: 51cm x 51cm x 10cm depth / 20" x 20" x 4" depth</b>
<b>Cushion Weight</b>	2 kgs

## FEATURES

Contoured cushion surface

Internal air channels

Waterproof, multi stretch, vapour permeable PU cover, impregnated with Acti-fresh biocide

Cover has white inner substrate

Cover can be removed and laundered at 71°C

Easy cleaning with neutral detergent and a single use wipe

All parts replaceable

Available with black cover

## BENEFITS

Provides greater weight distribution and effectively reduces pressure

Provides increased air circulation and reduces the build up of heat and moisture

Reduces shear and friction forces whilst aiding infection control

Helps to aid cushion audits by making staining easily identifiable

Aids decontamination protocols

For ease of cleaning if full decontamination is not required

Provides cost effectiveness

Suitable for a wide range of users in differing environments

\*All equipment should be used in conjunction with clinical judgement and assessments dependent on individual patient morphology. Essential nursing care is pivotal in providing the optimum pressure ulcer prevention plan.