


Bariatric Bed

- Full-electric, three-function operation
- New DC Electronics that are quieter and offer a faster Hi-Lo function than traditional AC
- DC Electronics are splashproof and offer a lower voltage making them more reliable and less prone to electronic failures
- Designed to meet the requirements of applicable HCPCS codes
- Split-frame with removable ends for easy transport
- Heavy-duty grid sleep deck
- Optional Swing-down style trigger release rail requires no tool assembly and designed with quick release for easy removal
- High-strength, steel-reinforced, truss frame
- Durable powder-coated finish
- Simple assembly: no special tools required
- Self-wrenching nut: assembled for easy part replacement
- Ball-screw bed ends
- Bed height: low 17" to high 27"
- Maximum patient weight is 600 lbs
- 700 lb evenly distributed total weight capacity including the patient, mattress, rails and accessories
- Overall bed dimensions: 46.5" width (with rails) and 87" length
- Sleeping surface dimensions: 42" width and 80" length
- Total Bed Weight: 295 lbs
- Also available as a package with 650-1501-1633 bariatric mattress and A8340 head half rails

Warranty:

- Two year warranty for defects in workmanship and materials of mechanical components, frame, and electronics
- During the warranty period, defective items will be repaired or replaced at manufacturer's option at no charge

ABL-B700	Full-Electric Bed Only	1 ea
ABL-B700-PKG	Bed, Head Half Rails, Bariatric Mattress	1 ea
650-1501-1633	Bariatric Mattress	1 ea
A8340	Head Half Rails	1 ea
A8510	Foot Half Rails	1 ea



5" locking caster