

STERLING

mobility



ELITE XS

Comfort, performance and style

The Sterling Elite XS has been designed for those who want to go just that little bit further.

An 8 mph scooter with all-round suspension and a fully adjustable seating solution, the Elite XS offers the ultimate in ride comfort, performance and style.

Designed to meet the needs of the most demanding user, the outstanding performance and comfort of the Elite XS offers you the ability to enjoy an active lifestyle.



ELITE XS



Comfort


- Fully adjustable sliding seat with extra lateral support
- Infinitely adjustable tiller
- Wigwag, finger or thumb control

Performance

- Top speed of 8 mph (12 km/h)
- Range up to 36 m (57 km)
- Adjustable acceleration and deceleration
- Rapid charging

Style

- Modern, sleek design
- Sporty alloy wheels
- High lustre paint finish



Control Panel Features

- Easy to read display
- Soft-touch buttons
- Horn
- Indicators
- Hazard warning lights
- 10 speed positions
- Battery level indicator

Steering Column Features

- Key for extra security
- Infinite tiller adjustment
- Parking brake
- Wigwag control
- Indicators
- Large removable basket
- Reflector
- Charging plug

Other Features

- Lockable storage compartment
- Powerful 18 W headlights
- Side and rear reflectors
- Easy to clean floor mat
- Rear lights
- Freewheel mechanism with audio signal
- Electromagnetic brake
- Anti-tip wheels
- Maintenance-free batteries
- Scratch resistant finish
- Brake light

Technical specifications

	3 wheel	4 wheel
maximum range*	57 km	57 km
forward speed maximum	12 km/h	12 km/h
reverse speed maximum	5 km/h	5 km/h
turning circle	115 cm	162 cm
incline 1:4 (25%), max. user weight	115 kg	115 kg
incline 1:5 (20%), max. user weight	150 kg	150 kg
total weight (excl. batteries)	82 kg	87 kg
total length	138 cm	138 cm
total width	67 cm	67 cm
weight heaviest part	64 kg	69 kg
batteries (max.)	74 Ah	74 Ah
weight batteries each (74 Ah)	26.5 kg	26.5 kg
adjustable seatheight	48 - 55 cm	48 - 55 cm
backrest height	52 cm	52 cm
seat depth	43 cm	43 cm
seat width	51 cm	51 cm

* Range is calculated under test conditions. Weight, terrain and weather conditions can affect range. There could be options and features depicted in this brochure that may not be available in your country. Please refer to your order form or contact your local Sunrise Medical dealer for further information.



Sunrise Medical Limited
High Street Wollaston
West Midlands DY8 4PS
England
Phone +44 (0) 1384 446622
Fax +44 (0) 1384 446644
www.sunrisemedical.co.uk



Distributor: