

Dual Release Folding Walkers

- Durable 1" aluminum tubing provides maximum strength while remaining lightweight.
- Raised H-frame design facilitates over-toilet positioning.
- Contoured, textured handgrips offer a safe, secure hand hold.
- Dual release folding mechanism enables user to fold walker sides independently.
- Adjustable height fits a broad range of user heights.
- Accommodates all leg attachments.
- Folded depth 4".
- Designed to meet the requirements of HCPCS code: E0135.
- Includes laminated color operating instructions.
- Limited lifetime warranty.

616070A Adult 1 ea
616170A Junior 1 ea



Item	Width Inside Handgrips	Width Inside Rear Legs	Overall Width	Depth Front to Rear	Minimum Height	Maximum Height	Approx User Height	Weight Capacity
616070A	16.5"	20.125"	22.75"	16.5"	32.5"	41.5"	5'4" - 6'2"	300 lb
616170A	16.5"	19.5"	22.5"	15.5"	27.5"	34.5"	4'7" - 5'3"	300 lb

