

LET YOUR IMAGINATION GO AND THE QUICKIE 2 WILL GET YOU THERE.



QUICKIE 2

Whether you're going down the street or across the continent, the Quickie 2 is the kind of chair you can take with you at a moment's notice. Fold it up, and it's flat enough to travel in places most other chairs can't. It's lightweight and easy to handle. And when you lock out the frame, it gives a rigid-like performance.

And if that isn't enough, the Quickie 2 is now available with a rounded rear frame and contoured axle plate. Add that to the trillion other different configurations available on the Quickie 2 and what you get is the most versatile chair on the market.

So Get Out There with the Quickie 2. The chair that works - and plays - as hard as you do.

QUICKIE
The RealChair People.



QUICKIE

2



- Weight: 27 lbs.
- Color Options: Blue, Black, Red, Midnight Purple, Silver, Sky White, Teal, Hot Pink, Ultra Yellow, Lavender, Blue Sapphire, Black Diamond, Candy Red
- Lifetime frame warranty



Depth Adjustable Back

The depth adjustable back provides 18° angle adjustment and is placed 2 1/2" farther back than a standard back, with forward adjustment of 1". The back works well in conjunction with alternative seating systems.



One-Arm Drive

The one-arm drive unit can be mounted on the chair's right or left side, is easily folded and adds less than 5 pounds to the weight of the chair.



Hanger Options

Quickie 2 is available with a variety of hanger options. Choose from 60°, 70° and 90° standard hangers, or try our 70° "V" style hanger.

QUICKIE 2HP

- Weight:
Quickie 2HP: 22 lbs.
Quickie 2HP Titanium: 20 1/4 lbs.
- Color Options: Blue, Black, Red, Midnight Purple, Silver, Sky White, Teal, Hot Pink, Ultra Yellow, Lavender, Blue Sapphire, Black Diamond, Candy Red, Black Splatter, Red Splatter, Sky White Splatter
- Titanium frame option
- Lifetime frame warranty



Platform Flip-Up Footplate

This one-piece, swing-away footrest locks into place and latches under the footrest tube when down. The footrest is also angle-adjustable to accommodate different front frame angles.



Performance Wheel Options

Choose from a variety of tires, rims and hubs to create the ultimate high performance wheels.



Rounded Rear Frame and Contoured Axle Plate

Contoured axle plates are available on the Quickie 2 and Quickie 2HP rounded rear frame option. This new design allows a 1 1/8" lower seat height.

QUICKIE 2

A. FRONT SEAT HEIGHT (For standard frames only)

Caster Fork Mounting Hole	5" LP Poly 5 1/4" Fork	5" LP Poly 6" Fork	5" LP Poly 7" Fork	6" Pneu 5 1/4" Fork	6" Pneu 6" Fork	6" Pneu 7" Fork	6" Poly 5 1/4" Fork	6" Poly 6" Fork	6" Poly 7" Fork	8" Pneu 5 1/4" Fork	8" Pneu 6" Fork	8" Pneu 7" Fork
Top	16 15/16"	17 3/4"	18 3/4"	17 1/2"	18 5/16"	19 5/16"	17 3/8"	18 5/16"	19 3/16"	18 1/2"	19 5/16"	20 5/16"
Bottom	17 15/16"	18 3/4"	19 3/4"	18 1/2"	19 5/16"	20 5/16"	18 3/8"	19 5/16"	20 3/16"	19 1/2"	20 5/16"	21 5/16"


B. REAR SEAT HEIGHT (For standard frames and axle position only)

Rear Axle Hole Positioning	20" Pneu	22" FP Poly	22" Pneu	24" LP Poly	24" FP Poly	24" HP Clincher	24" Pneu or Kevlar	26" HP Clincher	26" Pneu
1&2	18 7/8"	19 1/2"	19 3/4"	20"	20 1/4"	20 3/8"	20 1/2"	21 1/4"	21 5/8"
2&3	17 3/4"	18 3/8"	18 5/8"	18 7/8"	19 1/8"	19 1/4"	19 3/8"	20 1/8"	20 1/2"
3&4	16 5/8"	17 1/4"	17 1/2"	17 3/4"	18"	18 1/8"	18 1/4"	19"	19 5/8"
4 & 5	15 1/2"	16 1/8"	16 3/8"	16 5/8"	16 7/8"	17"	17 1/8"	17 7/8"	18 1/4"

C. OVERALL WIDTH

Frame Width	12"	13"	14"	15"	16"	17"	18"	19"	20"	Folded
Mag wheel with aluminum handrims	19"	20"	21"	22"	23"	24"	25"	26"	27"	12 1/2"
Spoke wheel with aluminum handrims	20"	21"	22"	23"	24"	25"	26"	27"	28"	12 1/2"

Actual measurements may vary up to 3/8". Front seat height must be equal to or greater than rear seat height. Front seat height cannot be more than 2" higher than rear seat height. Please consult your order form or customer service for combinations and restrictions. Specifications subject to change without notice.



QUICKIE 2HP/2HP TITANIUM

A. FRONT SEAT HEIGHT (For standard frames only)

Caster Fork Mounting Hole	5" Micro 5 1/4" Fork	5" Micro 6" Fork	5" Micro 7" Fork	5" LP Poly 5 1/4" Fork	5" LP Poly 6" Fork	5" LP Poly 7" Fork	6" Pneu 5 1/4" Fork	6" Pneu 6" Fork	6" Pneu 7" Fork	6" Poly 5 1/4" Fork	6" Poly 6" Fork	6" Poly 7" Fork
Top	16 5/16"	17 1/8"	18 1/8"	17 3/16"	18"	19"	17 3/4"	18 9/16"	19 9/16"	17 5/8"	18 7/16"	19 7/16"
Bottom	17 5/16"	18 1/8"	19 1/8"	18 3/16"	19"	20"	18 3/4"	19 9/16"	20 9/16"	18 5/8"	19 7/16"	20 7/16"

B. REAR SEAT HEIGHT (For standard frames and axle position only)

Rear Axle Hole Positioning	22" Pneu	24" Lates Tubular	24" HP Clincher	24" Pneu	24" HT Tennis/B-ball Turbo Tread	26" Lates Tubular	26" HP Clincher	26" Pneu	26" HT Tennis/B-ball Turbo Tread
1&2	19 3/4"	20 1/4"	20 3/8"	20 1/2"	20 5/8"	21 1/4"	21 1/4"	21 5/8"	22"
2&3	18 5/8"	19 1/8"	19 1/4"	19 3/8"	19 1/2"	20 1/8"	20 1/8"	20 1/2"	20 7/8"
3&4	17 1/2"	18"	18 1/8"	18 1/4"	18 3/8"	19"	19"	19 3/8"	19 3/4"
4 & 5	16 3/8"	16 7/8"	17"	17 1/8"	17 1/4"	17 7/8"	17 7/8"	18 1/4"	18 5/8"

C. OVERALL WIDTH

Frame Width	12"	13"	14"	15"	16"	17"	18"	19"	20"	Folded
Mag wheel with aluminum handrims	19"	20"	21"	22"	23"	24"	25"	26"	27"	12 1/2"
Spoke wheel with aluminum handrims	20"	21"	22"	23"	24"	25"	26"	27"	28"	12 1/2"

Actual measurements may vary up to 3/8". Front seat height must be equal to or greater than rear seat height. Front seat height cannot be more than 2" higher than rear seat height. Please consult your order form or customer service for combinations and restrictions. Specifications subject to change without notice.



SUNRISE

Sunrise Medical Limited
High Street Wollaston
West Midlands DY8 4PS
England
Tel +44 (0) 1384 446688
Fax +44 (0) 1384 446699
www.sunrise-medical.co.uk

