

Neon

First class travel...



Neon

The **Mini-Max** Principle

Compact, folding and lightweight for easy stowage in a plane, train or car.

Lightweight frame weighs just 7.1 kilos with intuitive folding cross brace with retaining strap, removable legrests and fold down backrest.



Neon Fixed Front version



Neon swing-away version

Getting free and independent



First class travel through fine adjustment... journey after journey.

Infinite fine tuning possibilities to ensure perfect driving ergonomics can be achieved and maintained.

Fashionable style and looks in whatever combination you wish to travel.

Up to 3 million option combinations with over 20 design choices.

3.7 million colour combinations including 6 upholstery colours and 27 brilliant and matt paint finishes for your frame, hubs, forks, wheels, handrims and footplates.



Modern open frame design (Fixed front)



Auto-folding footboard (option)



Small folding (Fixed front)



Internal cross-brace ensures optimum rigidity and drive performance



Simple seat height adjustment with preset factory "centre of gravity" settings



Improved, durable folding backrest system offering angle adjustment from -12° to $+8^{\circ}$



Neon Swing-away: with swing-in and swing-out function




Fast folding frame and back with quick release legrest ensures convenient transportation



Wide range of side guards and armrest options

There is room for the **Neon**
even in the smallest car.



Technical data sheet

27 individual frame colours		Folded width:	37cm (14.5") (SA) / 38cm (15") (FF)
6 individual sling colours		Optional seat depth growth	2.5 and 5cm (1 and 2")
Max. user weight:	115kg* (18 stone) / 140kg** (22 stone)	Weight of the heaviest part:	7.4kg (SA) / 7.5kg (FF)
Total weight:	From 9.9kg	Transit tested to ISO 7176 Part 19 and passed	Yes

* Seat width 42cm (16.5") ** Above seat width 42cm (16.5")

For further information on the full specification, options and accessories, please refer to the order form.

All information is subject to change without notification. Please consult Sunrise Medical with any queries you may have.

If you are visually impaired, this document can be viewed in PDF format at www.SunriseMedical.com or alternatively is available on request in large text.



PB_NEON_EN_GB_REV1

Sunrise Medical Limited
High Street, Wollaston
Stourbridge
West Midlands
DY8 4PS

Tel: +44 (0) 1384 44 66 66

Fax: +44 (0) 1384 44 66 44

www.SunriseMedical.com