



# GT/GTi Series

Innovative mono-frame chairs offering maximum strength and performance at the lowest possible weight.



# QUICKIE® GT

The flagship ultra lightweight aluminum chair

- Weighs in at less than 21.5 lbs



Quickie GT  
with XTR suspension



## Lifestyle Options

### Give your chair some attitude

- Over 20 paint colors to choose, including 8 new Kolorfusion™ colors
- Anodized parts
- New back release lanyard.

### Without sacrificing comfort.

- Natural Fit™ Handrims
- New Folding Lockdown Angle Adjustable Backrest
- GT with XTR suspension




Quickie GTi with KolorFusion Digital Camo

Natural-Fit Handrim



Catch-A-Wave



SpotFires



Mossy  
Oak Camo



Woodland  
Camo



Digital  
Camo



Digital  
Woodland  
Camo



Desert  
Camo



Digital  
Desert  
Camo





**You drive our design.** In development, we consider how you interact with the chair and how both can work as one to maintain maximum efficiency, performance, style and fit. Before you receive the chair, the product undergoes rigorous testing in addition to the required testing necessary for federal guidelines. Through our focus on these core elements: design, development and testing; we design these chairs for you and by you.

# QUICKIE® GTI



## Lightweight

Feel the responsiveness of a true ultra light chair

- Reactive to your movements
- Ease of propulsion
- Minimize Repetitive Strain Injury (RSI)
- Improved transportability
- Includes the LITE Spoke wheels as a standard option

## Adjustability

Fit the chair to you with customized adjustments

- Multiple caster and axle adjustments allows you to fit the chair to your environment
- New 3-D Perpendicular Caster Alignment technology.
- Designed for comfort through back recline correction and new Folding Lockdown Angle Adjustable Backrest

**The strongest, lightest titanium rigid chair in today's market.**

- Elite performance in an ultra lightweight design





## COMPLEMENT YOUR QUICKIE


Order your Jay® Cushion and Jay® J3 Back shipped complete on-chair with any Quickie GT/GTi chair. Just check the boxes on the order form and save time and money!




The engineers at JAY take the work of designing JAY cushions very seriously. Over the past 8 years they have commissioned 5 major scientific studies and continue to actively work with a number of top researchers in the academic community. To the untrained eye it may look like a simple cushion, but to our engineers it is a synthesis of ideas from the best minds in the world.

Here at JAY, we make the **Science of Seating** our motto.

### Jay Cushions



**Basic**  
E2601/E2602



**Soft Combi P**  
E2605/E2606



**Easy**  
E2607



**Lite**  
E2603/4  
E2607/8



**J2**  
K0734/5  
K0736/7



**J2 DC**  
K0734/5  
K0736/7

### Jay Backs



**J3 Posterior**  
E2613



**J3 Posterior & Lateral**  
E2615



**J3 Posterior & Deep Lateral**  
E2620

## Quickie GT/GTi Series Specifications

SPECIFICATIONS	GT	GTi
Frame Style	Rigid – open	Rigid – open
Transport Weight	12.4 lbs.	10.9 lbs.
Overall Weight	21.3 lbs.	17.9 lbs.
User Weight Capacity	265 lbs.	265 lbs.
Seat Widths	12" -20"	12" -20"
Seat Depth	12" -22"	12" - 22"
Back Angle Adjustment	85° – 102°	85° – 102°
Front Seat To Floor Height	16" - 21"	16" - 21"
Rear Seat to Floor Height	13" - 20"	13" - 20"
Frame Angle	75°, 80°	75°, 80°, 85°
<b>WARRANTY</b>		
Warranty	Limited/Lifetime	Limited Lifetime



Customer Service: 800.333.4000

For additional product specifications please visit:  
[www.sunrisemedical.com](http://www.sunrisemedical.com)

©Sunrise Medical, Inc. 02.09  
480032 Rev. A