

QUICKIE®

Easily and fully adjusts to your child's growing requirements

**15% more power
for 30 km* range!**
Easier servicing!

Samba²
Kid



The Samba² Kid is fun and active with powerful batteries and centre of gravity adjustments, ensuring the best autonomy and driving response. Adjustable seat width, depth, length and height ensures the chair grows with your child. Manual seat and back angle adjustment ensures comfort, and the Samba² Kid de-docks and disassembles without tools for easy travel.

* 30 km range when 60 Ampere batteries chosen as option (ISO 7176-4)



Dimensions:



Seat height and centre of gravity adjustment (one tool).



Kerb Climber (option) and 9 inch castors (standard).



Adjustable suspension for improved comfort and ride.




Hand controls to adjust seat and back angle.



Speciality Controls to meet the full range of paediatric requirements.



De-docks and disassembles without tools for easy transport.



Range 60 Ah Batteries: ISO 7176 - 4	30 km
Speed	6 km/h (10 km/h optional)
Kerb climbing	max 5 cm / max 10 cm with kerb climber option

Turning radius	80 cm
Battery	2 x 50 Ah (60 Ah optional)
Max safe slope	10° (18%)

STEPS™
Sunrise Training and Education Programmes
KNOWLEDGE FOR THE FUTURE

SUNRISE MEDICAL®

Sunrise Medical Limited
High Street
Wollaston, Stourbridge
West Midlands DY8 4PS, England
Phone +44 (0) 1384 446622
Fax +44 (0) 1384 446644
www.sunrisemedical.com