



**Youngster 3**



**QUICKIE**

**INNOVATE.  
BUILD FOR THE FUTURE**

**Discover new places,  
make new friends, have fun  
and play!**



**Youngster 3 –  
the new lightweight with  
superb design!**

# QUICKIE

The latest technology combined with proven solutions guarantee that the new Youngster 3 meets the needs of parents, therapists and users alike.

## Lightweight by design

The lightest overall weight with optimum functionality was the most important aim when the Youngster 3 was being developed. The result is a folding children's wheelchair which will impress you and, above all, in comparison to the previous version gives enormous weight-savings on the castor mounting, rear wheel axle and the cross-brace.

Even with the added options, the smallest amount of weight has been saved. So the Youngster 3 is available with lightweight rear wheels, with breathable seat sling or with the carbon-fibre look sideguard.

## Superb design

A wheelchair is not only a living aid, but also an expression of individuality and lifestyle. Therefore, when the Youngster 3 was designed, particular attention has been paid to shape and design.



As well as a modern range of colours, the Youngster 3 stands out, above all, by its well-designed look. The complete front frame has been re-designed and this, together with the stylised forks and the new spokeguard, give the Youngster 3 an overall feeling of lightness.

## Functionality for maximum individuality

The new Youngster 3 can be quickly and easily adapted to suit individual requirements and is available in three frame variations, with different seat widths and depths.

## Options

The Youngster 2 always stood out from the crowd because of the wide range of options available. Building on these high standards, the Youngster 3 offers options such as:

- Transit Brackets
- Jay seat belt
- Swingaway mechanism for the platform and two-part foot board
- Lightweight options
- Seat shell rear frame depth extension (5 cm)

# QUICKIE

## Youngster 3



### Dimensions

■ Seat Widths	22 cm to 40 cm (in 2 cm steps)
■ Seat Depths	26 cm to 40 cm (in 2 cm steps)
■ Seat Heights	38 cm to 49 cm
■ Lower Leg Lengths	10 cm to 45 cm (depending on configuration)
■ Back Heights	25 cm, 30 cm, 35 cm, 40 cm (+/-2.5 cm)
■ Growth	enabled via growth kits

#### ■ Weight

- Frame with platform style footrest 7.7 kg
- Aluminium sideguard with clothes protector 0.7 kg
- Carbon-fibre look sideguard with integrated clothes protector 0.4 kg
- Standard rear wheels with aluminium handrims 3.2 kg
- Lightweight rear wheels with aluminium handrims 2.7 kg
- Standard rear wheels with stainless steel handrims 3.8 kg

Rear Wheels	22" and 24"
Rear Wheel camber	1, 2, 3, 4 degree option
Castors	4", 5" and 6"
Maximum user weight	75 kg



Medical Device Directory applied for: Youngster 3 Item No.: 18.50.03xxxx

**ST<sub>PS</sub>**  
Sunrise Training and Education Programmes

**KNOWLEDGE FOR THE FUTURE**

 **SUNRISE  
MEDICAL**



Sunrise Medical Limited  
High Street Wollaston  
West Midlands DY8 4PS  
England  
Phone +44 (0) 1384 446622  
Fax +44 (0) 1384 446644  
[www.sunrisemedical.com](http://www.sunrisemedical.com)

**INNOVATE.  
BUILD FOR THE FUTURE**