Zippie Colors



Stars and Stripes

Zebra* * Black Zebra Kolorfusion™ pattern is available on any base color. Shown here on Candy Blue base color.

NOTES:

Print processes cannot accurately represent paint finishes.

Actual colors may vary.

Colors and patterns may vary based on type of chair or frame material.

Not all colors available on all products. Please see product order form for specific color options.



Zippie TS - Shown with NEW "Glow" paint.



Customer Service: 800.333.4000 www.SunriseMedical.com

©2010 Sunrise Medical, Inc. 05.10 103731 Rev. C





Pediatric Wheelchair Reference Guide

zippie®	FOLDING W/TILT	FOLDING		RIGID W/TILT		RIGID		
	Zippie TS / TS SE	Zippie GS / GS SE	Zippie 2	Zippie TS / TS SE	Zippie IRIS / IRIS SE	Zippie ZONE	Zippie Kidz	Zippie GS / GS SE
			OFTER					
Key Selling Points	 Tilt-in-space up to 45° Lightweight welded aluminum frame construction ideal for portability Designed to accept the widest possible variety of seating solutions 	 Built-in growth system Available in reverse configuration Proven design and the highest quality construction for trouble-free durability 	 Locking cross brace gives it an easy-pushing solid feel Extremely adjustable with the ability to grow as the child grows Swing-away footrests and quick-release components One-Arm Drive option 	 Tilt-in-space up to 45° Lightweight welded aluminum frame construction ideal for portability Designed to accept the widest possible variety of seating solutions 	 Patented rotation-in-space technology ensure center of gravity always remains in the center of the base Lightweight Can be disassembled for transport 	 Ultra lightweight Growth capabilities Easy adjustments Cool Style 	 Large wheels can be placed in front or rear Width and depth adjustments allow the chair to grow with the rider Easy chair adjustments and set-up 	 Built-in growth system Available in standard configuration (Reverse configuration shown) Proven design and the highest quality construction for trouble-free durability
HCPCS CODE	E1234	E1238	E1238	E1233	E1233	E1235	E1235	E1237
Transit Option	Yes	Yes	Yes	Yes	Yes	No	No	Yes
Tilt	45°			45°	-5° to 50°, 0° to 55°, +5° to 60°, -5° to 35°, 0° to 40°, +5° to 45°			
Chair Weight*	29 lbs w/o footrests	25 lbs. w/o footrests	25 lbs. w/o footrests	29 lbs. w/o footrests	22 lbs.	14 lbs.	14 lbs. w/o footrest	25 lbs. w/o footrest
Weight Capacity (w/o seating)	165 lbs.	165 lbs.	165 lbs.	165 lbs.	I 65 lbs.	165 lbs.	I 50 lbs.	165 lbs.
Overall Width*	18.25" to 24.25"	19" to 27"	20" to 24"	18.25" to 24.25"	Frame width +8"	Seat width +7" •	18" to 22"	19" to 27"
Tubing Diameter	Frame: I" Seat Rail: I"	Frame: I" Seat Rail: I"	Frame: I" Seat Rail: I"	Frame: I" Seat Rail: I"	Frame: I" Seat Rail: I"	Oval	Frame: I" Seat Rail: I"	Frame: I" Seat Rail: I"
Frame Seat Width	10" to 18"	10" to 18"	2" to 6"	10" to 18"	10" to 18"	8" to 16"	10" to 14"	10" to 18"
Seat Depth	I3" to 20"	10" to 20"	2" to 8"	3" to 20"	10" to 18"	8" to 18"	8" to 14"	10" to 20"
Back Angle Adjustment	85° to 120°	85° to 120°	85° to 120°	85° to 120°	85° to 120°	–8° to 27°	90°, 96°, 102°	85° to 120°
Front Seat-To-Floor Height	15.5" to 20"	14" to 20"	15" to 20"	15.5" to 20"	12.5" to 19.5" w/drop seat: 11" to 17.5"	14" to 19"	10 3/4" to 14"	14" to 20"
Rear Seat-To-Floor Height	4" to 9"	12.5" to 20"	13.5" to 19"	4" to 9"	12.5" to 19.5" w/drop seat: 11" to 17.5"	13" to 19"	9 3/4" to 14 1/4"	12.5" to 20"
Footrest / Frame Angle	60°, 70°, 80°, Swing-In/Out, 70°, 80°, 90° Swing-Away, HD, ELR, ALR	60°, 70°, 80°, Swing-In/Out, 70°, 80°, 90° Swing-Away, HD, ELR, ALR	60°, 70°, 80°, Swing-In/Out, 70°, 80°, 90° Swing-Away, HD, ELR, ALR	60°, 70°, 80°, Swing-In/Out, 70°, 80°, 90° Swing-Away, HD, ELR, ALR	60°, 70°, 80°, Swing-In/Out, 70°, 80°, 90° Swing-Away, HD, ELR, ALR	80°	90°	60°, 70°, 80°, Swing-In/Out, 70°, 80°, 90° Swing-Away, HD, ELR, ALR
Back Type	Fixed Height and Height Adjustable	Fixed Height and Height Adjustable	Fixed Height and Height Adjustable	Fixed Height and Height Adjustable	Fixed Height and Height Adjustable	Folding Angle Adjustable with double lock-down	Folding, Rigid	Fixed Height and Height Adjustable
Back Height	Fixed Stroller Handles 19", 22", 25", Height Adjustable (16"to22", 19"to25")	Fixed Stroller Handles 19", 22", 25", Height Adjustable (16"to22", 19"to25")	Fixed Stroller Handles 19", 22", 25", Height Adjustable (16"to22", 19"to25")	Fixed Stroller Handles 19", 22", 25", Height Adjustable (16"to22", 19"to25")	Fixed Stroller Handles 19", 22", 25", Height Adjustable (16"to22", 19"to25")	6" to 8" up to 14" to 18"	Regular (8"to10"), Tall (10"to12"), Stroller Handle	Fixed Stroller Handles 19", 22", 25", Height Adjustable (16"to22", 19"to25")
Center of Gravity Adjust. (Varies Depending on Configuration)	No	Horizontal – Yes Height – Yes	Horizontal – Yes Height – Yes	No	Horizontal – Yes Height – Yes	-0.25" to 2.25"	No	Horizontal – Yes Height – Yes
Rear Wheel Sizes	12", 16", 18", 20", 22", 24"	12", 16", 20", 22", 24"	12", 16", 18", 20", 22", 24"	12", 16", 18", 20", 22", 24"	12", 16", 20", 22", 24"	18", 20", 22", 24"	12", 16", 20"	12", 16", 18", 20", 22", 24"
Caster Sizes	3" to 8"	3" to 8"	3" to 8"	3" to 8"	4" to 8"	3" to 6"	3", 5"	3" to 8"
True Fit or One-Time Free Growth Kit	Z TS: Yes Z TS SE: No	Z GS: Yes Z GS SE: No	Yes - Growth	Z TS: Yes Z TS SE: No	Z IRIS: Yes Z IRIS SE: No	No	Yes - Growth	Z GS: Yes Z GS SE: No
Warranty**	Frame: Lifetime Parts: I Year	Frame: Lifetime Parts: Year	Frame: Lifetime Parts: Year	Frame: Lifetime Parts: 1 Year	Frame: Lifetime Parts: I Year	Frame: Lifetime Parts: 1 year	Frame: Lifetime Parts: 1 Year	Frame: Lifetime Parts: Year
Intended Consumer	Juvenile 4 – 16 years old. Tilt. Accepts seating. Transportability	Juvenile 4 – 16 years old. Self propeller.Accepts seating. Reverse wheel configuration opt.	Juvenile 4 – 16 years old. Folding. Accepts seating. Self propeller.	Juvenile 4 – 16 years old. Tilt. Accepts seating.	Juvenile 4 – 16 years old. Tilt (Rotation). Accepts seating.	Juvenile 4 to 16 years old. Self propeller.	Toddler I – 3 years old. Self propeller. Reverse wheel configuration opt.	Juvenile 4 – 16 years old. Self propeller. Accepts seating. Reverse wheel configuration opt.

Specifications are subject to change without notice.

* Actual measurements may vary according to the chair configuration and accessories.

** Refer to Owner's Manual for details.