

zippie®

TS™

Pediatric Folding Tilt and Positioning



Featured With New! Glow in the dark paint

The **Zippie TS** is the original premier pediatric tilt-in-space chair that is available in either a folding or rigid style frame. The Zippie TS has built-in frame growth with the widest variety of options and is able to accept many seating and positioning solutions.

Tilt-in-Space up to 45 degrees. Available in both rigid and folding frames. Also available with **Foot Release Tilt Actuator**, which eliminates exposed cables.

Ability to grow the chair frame as the child's needs change. The Zippie TS™ provides growth in both seat depth and seat width.

Optional Folding Frame is a lightweight welded frame construction which allows the frame to compress into a small package that is ideal for transportability and storage in tight places.



ZIPPIE TS SPECIFICATIONS

	Folding w/Tilt	Rigid w/Tilt
HCPCS Code:	E1234	E1233
Overall Weight*:	29 lbs. w/o footrest	29 lbs. w/o footrest
Weight Capacity:	165 lbs.	165 lbs.
Overall Product Width:	18.25" -24.25"	18.25" -24.25"
Transit Option:	Yes	Yes
Front Seat Height:	15.5" - 20"	15.5" - 20"
Rear Seat Height:	14" - 18.5"	14" - 18.5"
Hanger Options:	60° S-A, 70°, 80°, 90° (S-A, Front Mount & HD), ELR, ALR	60° S-A, 70°, 80°, 90° (S-A, Front Mount & HD), ELR, ALR
True Fit Growth System:	Yes	Yes
Warranty:	Lifetime on frame & cross brace, 1 year on parts	Lifetime on frame & cross brace, 1 year on parts

* Actual measurements may vary according to the wheelchair's configuration and accessories