

## QUICKIE IRIS GROWTH GUIDE (EIZ4)

When growing in **width**, the following parts of the chair are possibly affected:

- Seat Frame
- Seat Pan
- Drop Seat Pan (width only)
- Removable Stroller Handle Extension
- Footplates/Heel Loops
- Footboard/Leg Strap
- One-Arm Drive Shaft
- Tray Table
- Shaft Growth Handle (Height adjustable backrest)
- Contracture Footrest Rail
- If upgrading to Heavy Duty, you will need one more Cross-Strut Tube for Base Frame.
- Vent Tray Strut Tubes (2 required)

When growing in **depth**, the following parts of the chair are possibly affected:

- Armrests
- Seat Pan
- Footrest Extension Tubes (Contracture Tubes also)
- Seat Rail (short, medium, long)
- Push Handles (Height growth: applies to Height Adjustable, Reclining, and Stroller Handles)
- Tilt Handle/Cable Assembly (Only if back height changes)
- Attendant Wheel Lock Cable

