



Rumba®



Rumba®

Compact, transportable, and adjustable....

Very compact base for lifts and narrow doors

The Quickie Rumba base is only 59cm wide (42cm seat frame width), allowing easy access into lifts and narrow doorways. This frame size can be adjusted up to 48cm client width via the armrest. For clients with the need of a 48cm seat frame the Rumba is still only 65cm wide.



Adjustable and Detachable

Easy to transport

Easily detachable back and lift out batteries combined with a folding frame allow easy transport, with the heaviest component weighing only 37.4kg. Ideal for storage in car boots and confined spaces.



Folding frame, detachable back, lift-out batteries

Easy to adjust seat, backrest and armrests

Can be adapted to clients from 42cm to 54cm width. Frame has 6cm armrest width adjustment.

Tool-free backrest height adjustment to 43, 46 and 49cm and up to 30° backrest angle adjustment complete with adjustable upholstery (option).

Detachable armrests with height adjustment. 3° of seat angle adjustment



Adjustable width



Height and angle adjustable backrest

Available kits

- Attendant control
- Lever drive attendant control
- Dual control
- Tension adjustable upholstery (grey)
- Anatomic padded seat
- Anatomic padded back
- Angle/depth adjustable footplates
- Standard headrest
- Anatomic headrest
- Swing-over tray
- Amputee support
- Crutch holder




Kerb climber (standard)

Rumba®

Compact, transportable, and adjustable....

Technical specifications



Max. speed	6kph / 3.7mph	Max. safe slope:	10°
Batteries:	40Ah / [50Ah]	Reversing width:	100cm (39.3")
Colours:	Blue, red, black and silver	Max. height of kerb:	100cm (39.3")
Max. user weight:	125kg (19.7 stone)	Max range:	30km (18.6 miles)
Total weight (with batteries):	75kg	Folded width:	40cm (15.7")
Weight of heaviest part:	37.4kg (frame)	Folded length:	83cm (32.6")
Turning circle:	85cm (33.4")	Folded height:	105cm (41.3")

For further information on the full specification, options and accessories, please refer to the order form.

All information is subject to change without notification. Please consult Sunrise Medical with any queries you may have.

If you are visually impaired, this document can be viewed in PDF format at www.SunriseMedical.com or alternatively is available on request in large text.



PB_RUMBA_EN_GB_REV1

Sunrise Medical Limited
High Street, Wollaston
Stourbridge
West Midlands
DY8 4PS

Tel: +44 (0) 1384 44 66 66
Fax: +44 (0) 1384 44 66 44

www.SunriseMedical.com