

Samba Updates

The Quickie Samba has been designed to suit the active person, right down to the last detail. As it is so easy to disassemble for effortless transportation, has a stable frame for excellent manoeuvrability and is completely equipped for the outdoors, the Samba meets all demands which may be required of a modern power chair, in day-to-day living and leisure time.

9" Front castors
(option)



Standard features

- Choice of 3 frame colours
- 70° Hangers
- Choice of left or right hand control mounting
- Composite footplates with heel loops
- VSI Integral control (Pilot+ as option)
- 12" Pneumatic rear wheels
- 40Ah Maintenance free batteries (50Ah as option)
- Seat widths from 12" to 20" available
- Base model available for JAY or other seating systems
- Manual tilt seat
- 220V/8A Automatic charger
- Manual recline back
- Transportation clamps
- Moulded height/width adjustable armrest (Detachable)
- Maximum user weight limit 120kg (19st)

Increased top speed
10kph*
(option)

Available from January 2005.

Samba kids ^{new}



The new Quickie Samba Kids may be smaller, but it shares all of the features of its bigger brother. Designed to grow with the user, the Samba Kids seat width can be easily adjusted from 12" to 16". Simple to use, excellent manoeuvrability and ideal for outside use, the Samba Kids allows children to be as creative and independent as they can be.

12" seat width



All information subject to change without notification. Please contact your local ASM for further information.



Sunrise Medical Limited
High Street Wollaston
West Midlands DY8 4PS
England
Phone +44 (0) 1384 446622
Fax +44 (0) 1384 446644
www.sunrisemedical.com

INNOVATE.
BUILD FOR THE FUTURE