



Quickie Sports & Recreation

Each of us possesses a spirit which is revealed in our STRENGTH, our ENERGY, our PASSION, our ABILITY, our CHOICE, our LIBERTY and our STYLE. What makes us all alike is the presence of a spirit. What makes us all unique is how we reveal it. Quickie Sports and Recreation is dedicated to providing innovative products to help you reveal your spirit.



Quickie Sports & Recreation



	C H A I R S			
Model	Ti All Court	All Court	Ti Match Point	Match Point
Material	Titanium	Aluminum	Titanium	Aluminum
Product Weight	21.5 lbs	24 lbs	20 lbs	22.5 lbs
Width	12"-20"	12"-20"	12"-20"	12"-20"
Depth/Length	10"-20"	10"-20"	10"-18"	10"-18"
Height	16"-21"	16"-21"	16"-21"	16"-21"
Weight Capacity	250 lbs	250 lbs	250 lbs	250 lbs
Camber Options	16° or 20°	16° or 20°	16° or 20°	16° or 20°
Front Seat Height	16"-21"	16"-21"	16"-21"	16"-21"
Rear Seat Height	13"-21"	13"-21"	13"-21"	13"-21"
# of Color Options	23	23	23	23
Anti-Tip Style	Single (std), Double	Single (std), Double	Single	Single
Support	Anti-tip brace	Anti-tip brace	Anti-tip brace	Anti-tip brace
Chest Straps	2", 4"	2", 4"	2", 4"	2", 4"
Wheel Type	Spinergy, Heat Treated, Cross Spoke	Spinergy, Heat Treated, Cross Spoke	Spinergy, Heat Treated, Cross Spoke	Spinergy, Heat Treated, Cross Spoke
Wheel Size	24", 26"	24", 26"	24", 26"	24", 26"
Key Option	Wing w/ standup	Wing w/ standup	extra wide 'footprint'	extra wide 'footprint'
Backrest	Folding, Non-Folding	Folding, Non-Folding	Folding, Non-Folding	Folding, Non-Folding



	S H A R K S	
Model	Shark	Shark S
Product Weight	33 lbs	30 lbs
Seat Width	14.8", 16.4", 17.9"	14.8", 16.4", 17.9"
Depth/Length	40.2"	40.2"
Height	8.6"	8.6"
Frame Width	15", 16.5", 18"	14"-18"
Weight Capacity	265 lbs	265 lbs
Rear angle adjust	30°-90°, 7 Basic, 3 fine settings	30°-90°, 7 Basic, 3 fine settings
Crank Length/Height	6.89"	6.89"-8.46"
Gear Speeds	7", 27"	27"
Gears	Sram S7™, Sram™5.0	Sram X9™
Brake	Shimano-Tiagra™	Shimano-Ultegra™
Back	adjust, foldable	adjust, foldable
# of Color Options	26	26

Shark Color Selecton



Customer Service: 800.333.4000
www.SunriseMedical.com

©10.2010 Sunrise Medical (US) LLC
 931218 Rev. H