

Breezy 300





Breezy 300

Simply light...



Weight...

One chair, one order form, multiple applications. The Breezy 300 provides a wide choice of backrests, rear wheels, accessories and seating options conforming perfectly to the individual needs of each user.

Functionality...

Attractive aluminium and composite design ensures that the Breezy 300 is a stylish, reliable and adaptable wheelchair.



Breezy 305



Breezy 315



Breezy 321



Breezy 312



Breezy 310



Breezy 341

Modern look...

Swing in, swing out detachable legrest, flip-up and detachable armrest and folding backrest options, ensure that the Breezy 300 is ideal for client transfers, accessing confined places and compact transportation.

Clients can be positioned via a range of headrests, anatomic cushions & backrests, elevating legrests and height adjustable armrests.



Long armrest (option)



Hemiplegic one arm drive with double handrim (option)



Drum brake for attendant use



Breezy 300

Simply light...

Breezy 300




370-520 mm




SW + 190 mm




410 mm




1060 mm



510 mm



450 mm



950 mm

Technical data sheet

Frame colour:	Brilliant black, selenium grey, apple green	Reversing width:	1070 mm
Upholstery colour:	grey	Folded width:	340 mm
Approved user weight:	120 kg	Folded height:	950 mm
Max. weight:	13 kg	Folded length:	1060 mm
Weight of heaviest part:	12.3 Kg		

Measures bases on model Breezy 305

For further information on the full specification, options and accessories, please refer to the order form.

All information is subject to change without notification. Please consult Sunrise Medical with any queries you may have.

If you are visually impaired, this document can be viewed in PDF format at www.SunriseMedical.com or alternatively is available on request in large text.



Kahlbachring 2-4
69254 Malsch/HD
Tel.: +49 (0) 72 53 / 9 80-0
Fax: +49 (0) 72 53 / 9 80-111

www.SunriseMedical.com