

Steel Transport Chair

The Everest & Jennings Transport Chair is designed as an affordable and safe patient transfer device. Its compact size and lighter weight allow a companion to lift and store the chair with ease.

Features:

- Available in Silver Vein Hammertone finish
- Fold down back makes transport and storage easy
- Black, padded, Nylon upholstery is flame retardant
- 17" and 19" seat widths
- Carry straps on seat allow for ease of lifting and folding
- Footrests are industry compatible and offer easy tool free adjustment
- Full length padded arms provide user comfort
- Weighs 24 lb w/o front rigging
- Weight Capacity 250 lb
- Designed to meet the requirements of Medicare code E1038

Warranty:

1 year on frame and crossbraces



Item Number	Description	Seat	Seat Depth	Seat Height	Back Height	Overall Depth	Overall Height	Overall Width	Folded Width	Weight
EJ795-1	Steel - Silver Vein	17"	15 3/4"	19"	18"	37"	39"	24"	10"	24 lb
EJ796-1	Steel - Silver Vein	19"	15 3/4"	19"	18"	37"	39"	24"	10"	24 lb

*NOTE: Everest & Jennings maintains the policy of continual product improvement, therefore, we reserve the right to make changes without notice. All dimensions are ±1/4". *Without Front Rigging.*

