

UniX



Unix

**Build to last...
The Fleet chair...**



Fleet management...

- Durable and reliable
- Common and interchangeable parts and accessories across entire XFamily (BasiX, UniX™, RubiX, HeliX and RelaX²)

 XFamily

Build to last...

Robust reliable design with full guarantee


Ensuring enduring mobility both inside and outside the home.



Flip-up sideguard with adjustable armpad long and short.



Rear wheel with maintenance free solid tyres and quick release axle



Standard padded Nylon; Fixed push handles

The Fleet chair...

Elevating legrests, legrests, height adjustable armrests and anti-tippers are just a few of the options and parts that are interchangeable across the new XFamily (UniX™, BasiX, RubiX, Relax²)



Elevating legrest (option)



Height adjustable armrest (option)




Tray options available. Slide on (pictured) and swing away.



UniX

Build to last...
The Fleet chair...


UniX




410-500 mm




SW + 190 mm




420 mm




1090 mm



510mm



430mm



970 mm

Technical data sheet

Frame colour:	Silver white	Weight of heaviest part:	12.2 kg
Upholstery colour:	black	Turning radius:	850 mm
Approved user weight:	125 kg	Max. Safe slope (with antitips):	10°
Max. weight:	18 kg		

For further information on the full specification, options and accessories, please refer to the order form.
All information is subject to change without notification. Please consult Sunrise Medical with any queries you may have.
If you are visually impaired, this document can be viewed in PDF format at www.SunriseMedical.com or alternatively is available on request in large text.



Kahlbachring 2-4
69254 Malsch/HD
Tel.: +49 (0) 72 53 / 9 80-0
Fax: +49 (0) 72 53 / 9 80-111
www.SunriseMedical.com