

## Preferred Care® Recliner Series, Tilt-in-Space (Model 565TG)

The Tilt-in-Space feature allows for improved resident positioning, increased comfort and pressure relief. Independent leg rest is easily activated by the caregiver from side of chair to provide leg elevation for increased resident comfort. Self storing, height adjustable drop arm is simple for the caregiver to activate and enables a safer, easier lateral transfer for the caregiver and resident. Padded armrests provide a soft surface for added resident comfort. Directional wheel lock on right rear caster reduces “crabbing” of recliner and enables the chair to be moved down hallways with less effort.

### Standard Features

- A unique suspension system self-adjusts to each person, cradling the resident from head to heel and virtually eliminating the risk of “bottoming out”.
- Extra-wide padded armrests and side panels maximize comfort.
- Durable epoxy coated welded steel frame is easy to clean and maintain.
- Quick-release upholstery allows for easy cleaning, removal and replacements. Vinyl upholstered back is sulfide stain, mildew and oil resistant.
- A variety of accessories are available to customize the chair for specific individual needs.
- Weight Bearing Capacity: 250 lb.

### • Lumex Wrap-Around Warranty

**Welded Frame:** Limited Lifetime

**Upholstery:** Limited One Year

**All Other Components:** Limited Three Year

<b>565TG427</b>	Blue Ridge	1 ea
<b>565TG851</b>	Doe Skin	1 ea
<b>565TG857</b>	Jade	1 ea
<b>565TG863</b>	Rosewood	1 ea

### (Meets California Technical Bulletin 133 Flammability Standards)

<b>FR565TG427</b>	Blue Ridge	1 ea
<b>FR565TG857</b>	Jade	1 ea
<b>FR565TG863</b>	Rosewood	1 ea



Blue Ridge (427)



Doe Skin (851)



Jade (857)



Rosewood (863)



- Upholstery color may vary slightly from brochure.
- Lumex Ortho-Biotic Seating with an FR prefix designates that the seating meets CA133 Fire Code specifications.
- Please contact Graham-Field Customer Service at 800-347-5678 for details on specialty colors and COM orders.

