

# Lift Chair Accessories

## Rise



Regal PR751TY in Night Navy  
Ultraleather® Brisa™ fabric shown.



Golden Technologies luxury lift chair line offers accessories to enhance the most comfortable lift chair made in the U.S.A. Make your luxury lift chair even more comfortable by choosing any of these fine lift chair accessories.

## Accessories



### Head Pillow

- Standard on PR-502
- Available as an option on all other models
- Style may vary



### Head Rest Cover

- Standard on all models



### Footrest Extension

- Extend the length of the chair by 5" when fully reclined
- Available on these models  
PR-401, PR-747, PR-752, PR-756MC, PR-505S, PR-505M, PR-505L, PR-501S, PR-501M, PR-501L, PR-501T, PR-501S-23, PR-501M-26, PR-501L-26, PR-501T-28, PR-355M, PR-355L, PR-359M, and PR-359L



## Additional Accessories

### ■ Armrest Covers

- Available on all models except the PR-410, PR-643, PR-751TY, PR-756MC, and PR-906
- Available in 13" or 18" lengths

### ■ Straight Lift Option

- This option lifts the user straight up out of the chair instead of at a forward angle.
- Available on the following models.  
PR-401, PR-747, PR-643, PR-410, PR-751TY, PR-752, PR-501JP, PR-501S, PR-501M, PR-501L, PR-501T, PR-501S-23\*, PR-501M-26D\*, PR-501L-26D\*, PR-501T-28D\*, PR-502\*, PR-200, PR-355M, PR-355L, PR-359M, and PR-359L

### ■ Left Hand Control Placement

- The control can be placed on the left side of the lift chair.
- Available on all models except the PR-751TY

### ■ Extra Pocket

- Available on all models

### ■ Far Infrared Heat:

Experience the relaxing effects of Far Infrared Heat, available on lift chairs exclusively from Golden! Far Infrared Heat is available on select models only. Consult your authorized Golden dealer or the Golden Customer Service Department for more information.

### ■ Deluxe Heat and Massage

- Enhance a lift chair with warming heat and a deep massage
- Available on all models except PR-200

### ■ Suede Fabric

### ■ Ultraleather® Brisa™ Fabric

### ■ Customer's Own Material

Customers can ship their own material to Golden Technologies for the perfect match to their decor! Lead times and upcharges will apply. Consult your authorized Golden dealer or the Golden Customer Service Department for more information.

### ■ Fabric Protectant

Select Golden fabrics come standard with a stain protectant included in the fabric. For the fabrics that do not come standard with stain protectant, Golden can apply that feature to the fabric of a lift chair before it is shipped from the factory. Consult your authorized Golden dealer or the Golden Customer Service Department for more information.

\*Straight lift option must be set at our factory prior to shipping for these models.



Golden Technologies • 401 Bridge Street, Old Forge PA 18518  
Tel: (800) 624-6374 • Fax: (800) 628-5165  
www.goldentech.com