

T R A D I T I O N A L S E R I E S


Rise



**Full chaise pad
now standard!**

Cambridge
PR-401

Softouch
with STAIN DEFENSE




Cognac




Sapphire



Sage



Mulberry



Rose



Williamsburg
PR-747



Pioneer
PR-643

Traditional Series

The Traditional Series features three models in classic styles that will subtly match any decor. The Cambridge and Williamsburg feature plush, overstuffed button backs for hours of comfort while the Pioneer offers a sleek, seam back and wood armrests for a style all its own.



Rise. Roll. Rest. Ride.



YOUR LIFE IN MOTION

TRADITIONAL SERIES



	Cambridge	Williamsburg	Pioneer
Model Number	PR-401	PR-747	PR-643
Size	Medium	Medium	Medium
Weight Capacity	375 Lb.	375 Lb.	375 Lb.
Recline Positions	3	3	3
Standard Fabrics	5	5	5
Fabric Options	28	32	32
Chaise Seat	YES	YES	YES
Back Type	Biscuit	Button	Seam
Head Rest Cover	YES	YES	YES
Arm Covers	YES	YES	NO
Battery Back Up	YES	YES	YES
Overall Width	33"	32.5"	29.5"
Overall Height	41.5"	44"	41"
Width Between Arms	20"	21"	21.5"
Floor to Top of Seat	19"	19.5"	19"
Seat Depth	20"	21.5"	21"
Seat to Top of Back	28"	28"	28"
Seat to Top of Arm	7"	6"	8"
Floor to Back Edge of Seat Raised	28"	28"	27.5"
Floor to Front Edge of Seat Raised	24"	24"	23"
Seat Angle Raised	20°	15°	20°
Floor to Top of Foot Rest Reclined	22"	23"	23"
Floor to Top of Back Reclined	33"	36"	34"
Back Angle Closed	110°	100°	115°
Back Angle Reclined	145°	130°	140°
Foot Rest Angle Closed	90°	90°	90°
Foot Rest Angle Opened	175°	180°	180°
Distance Required from Wall Reclined	20"	18"	18"
Extension reclined from Back of Seat to end of Foot Rest	41.5"	42"	40.5"
Life Time Warranty, Chair Frame, Lift Frame, Recline Mechanism	YES	YES	YES
Warranty on Electrical Parts	3 Yrs.	3 Yrs.	3 Yrs.
Pro Rated Years 4 Thru 7	YES	YES	YES
Warranty on Mechanical Labor	3 Yrs.	3 Yrs.	3 Yrs.
Warranty on Electrical Labor	1 Yr.	1 Yr.	1 Yr.
Removable Back Standard	YES	YES	YES
Smart Tek System Standard	YES	YES	YES
Customizable	NO	NO	NO
Left Hand Control Option Available	YES	YES	YES
Extra Pocket Option Available	YES	YES	YES
Heat & Massage Option Available	YES	YES	YES
Foot Rest Extension Option Available	YES	YES	NO
Straight Lift Option Available	YES	YES	YES

Sizes may vary depending on fabric, filling material, upholstery and carpet thickness. All measurements taken on concrete floor using levels and metal rulers.



Smart Tek

Our luxurious lift and recline chairs are equipped with electronics to help assure a secure lift every time you need it.



501 M shown

Golden Craftsmanship

Each luxury lift & recline chair is hand crafted in our facility using decades of motion furniture experience.



501 M shown

Added Comfort and Safety

Designed to perform year after year using advanced technology backed by the best warranty.



501 M shown

Safety System

Batteries not included. Safety System is designed to bring chair from recline to seated position in the event of a power outage.



Battery Backup System

Smart Engineering

Smooth, quiet lifting system offers a secure lift for the best performance.



Literature is current at time of printing. Golden Technologies reserves the right to make changes to the product or literature at any time.



Complies with UL 73 and CSA C22.2 No. 68

Listing # E112417

